



COMMERCIAL CENTER & R-D CENTER
F10, Building 3, 89 Shengli Road, Jiangning District 211100, Nanjing, China
Tel. 0086 25 86216803
Fax. 0086 25 86214463
E-mail. info@citotest.com

PRODUCTION
No. 339 Beihai West Road, Haimen, 226100 Jiangsu, P.R. China
Tel. 0086 513 82259348
Fax. 0086 513 82110195



Histology Pathology

CITOTEST®

Contents

MICROSCOPE SLIDES

- 6 Plain Microscope Slides
- 7 Frosted Microscope Slides
- 8 ColorCoat™ Microscope Slides
- 9-10 PrintAID™ Microscope Slides
- 11 Photosensitive Microscope Slide
- 12 ColorInx™ Microscope Slides adaptable to Etching Slide Printers
- 15 Color-PLUS™ PCB Microscope Slides
- 16-19 Color-PLUS™ PCL Microscope Slides
- 20-23 Color-PLUS™ PCI Microscope Slides
- 24 Liquid-Base Cytology Microscope Slides
- 26 Auto-Smear™ Microscope Slides
- 27 PapSmear™ Microscope Slides
- 27 Plastic Spatula For Blood Smear
- 28-37 Diagnostic Microscope Slides
- 38 Mega-Sized Microscope Slides
- 38 Microscope marking pen

MICROSCOPE COVER GLASS

- 40 EstaGlas™ Microscope Cover Glass
- 41 BoroGlas™ Microscope Cover Glass
- 42 EcoGlas™ Microscope Cover Glass
- 43 Circular Microscope Cover Glass
- 44 Hemacytometer Cover Glass
- 45 Non-Sticking Microscope Cover Glass

PROCESSING EMBEDDING CASSETTES

- 45-60 Processing/Embedding Cassettes
- 60-62 Base Molds/Biopsy Bags/Biopsy Foam Pads
- 62 Cassette Dispensing Sleeve Set
- 63 Paraffin Wax
- 63 Microtome Blades
- 64 Staining Rack/Staining Jar

TRANSPORT STORAGE ARCHIVING


- 66 Cardboard Tray, Mailer & Wallets
- 67-68 Plastic Tray, Mailer & Wallets
- 69-70 Slide Storage Box
- 71-73 Storage Cabinets
- 73 Slide Dispenser
- 74 Secure Archival Storage Solution

CITOTEST SCIENTIFIC CO.,LTD

CITOTEST was established in 1996 in East China. The owner ambitiously founded CITOTEST for creating a cost-effective Manufacturer and Provider for high quality Microscope Slides and Microscope Cover Glass. And based on a lot of fundamental hard work as well as persistent pursuit during the initial years, we finally made it and had CITOTEST well distinguished from other suppliers of economical products, and CITOTEST obviously became popular as a symbol of high quality products.

Being a private company, CITOTEST has been constructing its organizations exactly following the needs of developments since the day of its birth. All the CITOTEST members clearly understand that only continuous diligence and developments can assure CITOTEST as well as its members outstanding over the increasingly fierce competition. The young CITOTEST team, having gone through considerable challenges, is becoming mature every day. And with the adoption of more and more professional operational systems and tools, CITOTEST is becoming more and more professional for satisfying our customers.

Foreseeing the near future, the young CITOTEST with its growing team shall continuously develop our product portfolio to meet more comprehensive customer demands and commit ourselves to delivering our partners high quality products and exact solutions.



YOUR **FIRST** CHOICE

Exact solution to meet up with the demands in applications

Based on the dedicated research of the applications in industry and persistent teamwork of R&D and Engineering, any demand of a microscope slide for use in routine laboratory and clinical pathology can be satisfied with the CITOTEST® microscope slide family

Comprehensive Choices upon your quality consideration


With a great number of microscope slide specifications, CITOTEST® are providing different choices in quality to match up with your demands.



Different standards can be incorporated into production according to your requirements

The CITOTEST® microscope slides are manufactured according to the ISO8037/1 international standards.

- The ISO 8037/1 international standard
- BS7011:SECTION 2:1;1989
- Various customized demands and requirements



CSW (Super White Glass)

The CSW series Microscope Slides are made of pure white glass sheets that bring with ideal optical properties and stable molecular structure, insuring accurate experiments results and safe specimen preservation.

Each CSW series Microscope Slide is carefully controlled on line and strictly inspected before packaging and delivery.

Standard Class

The Standard Class Microscope Slides are produced of green flat glass sheets according to the CITOTEST® standards.

Super Grade

The Super Grade Microscope Slides are produced of 3/4 water white float glass sheets featuring high consistency in thickness and great surface evenness, delivering an added guarantee to instrument adaptability and experimental accuracy.

Each Super Grade serie Microscope Slide is carefully controlled on line and strictly inspected before packaging and delivery.

Choosing the CITOTEST® Super Grade Microscope Slides, you are able to experience ideal experiments with a cost-effective alternative apart from the CSW series Microscope Slides.

CITOGLAS.

The CSW series and Super Grade series Microscope slides

CITOPUS

The Standard Grade and Selected Microscope slides



STEPS FOR CHOOSING A CITOTEST MICROSCOPE SLIDE



A CUSTOMER DEMAND



APPLICATION

Define the intended applications of the slides that you need



EXPERIMENTAL EFFECTS

Consider the importance of your experimental results



COST

Make a cost-effective choice

THE RANGE

Regular Microscope Slides

ColorInx™ Microscope Slides

Adhesion Microscope Slides

Liquid-Based Cytology Microscope Slides

Auto-Smear™ Microscope Slides

PapSmear™ Microscope Slides

Plastic spatula for blood smear

Diagnostic Microscope Slides

Mega-Sized Microscope Slides

Microscope Slides with Cavity

Symbols


Printable by Thermo SlideMate SlidePrinters	Printable by Signature Slide Printers	Printable by Etching Slide Printers	Printable by Laser Slide Printers
Printable by Inkjet IPS Slide Printers	Adhesion microscope slides PCB series	Adhesion microscope slides PCL series	Adhesion microscope slides PCI series
Cut Edges 90°corners	Ground Edges 90°corners	Ground Edges 45°corners	Beveled Edges 45°corners

All the slides are precisely produced and controlled on our automatic production systems and the product quality consistency is highly ensured with the up-to-date engineering abilities as well as strict quality standards.



Regular Microscope Slides

Designed for routine experiments where no extra adhesion or adaptability to slide labelers is required in histology, pathology, hematology, cytology, microbiology and etc, for blood or body liquid tests, frozen tissue section examinations, biological specimen preparations, etc..



Plain Microscope Slides

Recommended for marking with the conventional labeling system, un-supportive for any direct printing or marking

Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

- Pre-cleaned to eliminate debris and improve wettability
- Beveled edges provide utmost safety during manual operations
- Super grade and super white glass meet higher transparency and more uniform thickness requirements




CATEGORY	CUT EDGES 90°CORNERS	GROUND EDGES 90°CORNERS	GROUND EDGES 45°CORNERS	BEVELED EDGES 45°CORNERS	PACKING
Super White Glass	0317-0001	0313-0001	0312-0001	0314-0001	50pcs/box, 1,000pcs/promo-pk,2,000pcs/carton
Super White Glass	0317-0021	0313-0021	0312-0021	0314-0021	72pcs/box, 1,440pcs/carton
Super Grade	0307-0002	0303-0002	0302-0002	0304-0002	50pcs/box, 1,000pcs/promo-pk,2,000pcs/carton
Super Grade	0307-0022	0303-0022	0302-0022	0304-0022	72pcs/box, 1,440pcs/carton
Standard Class	0307-0005	0303-0005	0302-0005	-	50pcs/box, 2,500pcs/carton
Standard Class	0307-0025	0303-0025	0302-0025	-	72pcs/box, 3,600pcs/carton

Frosted Microscope Slides

Recommended for marking with the conventional labeling system such as pencil or mark pens, Double frosted microscope slides create a considerable workload reduction over the Single Frosted Microscope Slides by offering a more random operation

Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

- Frosted area of approx. 20 mm depth, resistant to all common chemicals and routine stains that are utilized in laboratory
- Beveled edges provide utmost safety during manual operations
- Super grade and super white glass meet higher transparency and more uniform thickness requirements



Single Frosted Microscope Slides

CATEGORY	CUT EDGES 90° CORNERS	GROUND EDGES 90° CORNERS	GROUND EDGES 45° CORNERS	BEVELED EDGES 45° CORNERS	PACKING
Super White Glass	0317-2101	0313-2101	0312-2101	0314-2101	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
Super White Glass	0317-2121	0313-2121	0312-2121	0314-2121	72pcs/box, 1,440pcs/carton
Super Grade	0307-2102	0303-2102	0302-2102	0304-2102	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
Super Grade	0307-2122	0303-2122	0302-2122	0304-2122	72pcs/box, 1,440pcs/carton
Standard Class	0307-2105	0303-2105	0302-2105	-	50pcs/box, 2,500pcs/carton
Standard Class	0307-2125	0303-2125	0302-2125	-	72pcs/box, 3,600pcs/carton

Double Frosted Microscope Slides


CATEGORY	CUT EDGES 90° CORNERS	GROUND EDGES 90° CORNERS	GROUND EDGES 45° CORNERS	BEVELED EDGES 45° CORNERS	PACKING
Super White Glass	0317-2201	0313-2201	0312-2201	0314-2201	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
Super White Glass	0317-2221	0313-2221	0312-2221	0314-2221	72pcs/box, 1,440pcs/carton
Super Grade	0307-2202	0303-2202	0302-2202	0304-2202	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
Super Grade	0307-2222	0303-2222	0302-2222	0304-2222	72pcs/box, 1,440pcs/carton
Standard Class	0307-2205	0303-2205	0302-2205	-	50pcs/box, 2,500pcs/carton
Standard Class	0307-2225	0303-2225	0302-2225	-	72pcs/box, 3,600pcs/carton

ColorCoat™ Microscope Slides (Microscope Slides with 1 color frosted area at 1 surface End)

Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The microscope slides are designed with 1 silky color frosted area at 1 surface End; the color frosted area is 20mm in depth and resistant to all common chemicals and routine stains that are utilized in laboratory.

- * Recommended for marking with the conventional labeling system such as pencil or mark pens.
- * The microscope slides with beveled edges are highly recommended for insuring utmost safety during manual operations.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton

MARKING AREA	CUT EDGES, 90° CORNERS		GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS		BEVELED EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0317-6101-02	0307-6101-02	0313-6101-02	0303-6101-02	0312-6101-02	0302-6101-02	0314-6101-02	0304-6101-02
Green	0317-6101-08	0307-6101-08	0313-6101-08	0303-6101-08	0312-6101-08	0302-6101-08	0314-6101-08	0304-6101-08
Purple	0317-6101-09	0307-6101-09	0313-6101-09	0303-6101-09	0312-6101-09	0302-6101-09	0314-6101-09	0304-6101-09
Orange	0317-6101-10	0307-6101-10	0313-6101-10	0303-6101-10	0312-6101-10	0302-6101-10	0314-6101-10	0304-6101-10
Pink	0317-6101-12	0307-6101-12	0313-6101-12	0303-6101-12	0312-6101-12	0302-6101-12	0314-6101-12	0304-6101-12
Tan	0317-6101-15	0307-6101-15	0313-6101-15	0303-6101-15	0312-6101-15	0302-6101-15	0314-6101-15	0304-6101-15
White	0317-6101-16	0307-6101-16	0313-6101-16	0303-6101-16	0312-6101-16	0302-6101-16	0314-6101-16	0304-6101-16
Yellow	0317-6101-17	0307-6101-17	0313-6101-17	0303-6101-17	0312-6101-17	0302-6101-17	0314-6101-17	0304-6101-17

72/PK Packing: 72pcs/box, 1,440pcs/carton

MARKING AREA	CUT EDGES, 90° CORNERS		GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS		BEVELED EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0317-6121-02	0307-6121-02	0313-6121-02	0303-6121-02	0312-6121-02	0302-6121-02	0314-6121-02	0304-6121-02
Green	0317-6121-08	0307-6121-08	0313-6121-08	0303-6121-08	0312-6121-08	0302-6121-08	0314-6121-08	0304-6121-08
Purple	0317-6121-09	0307-6121-09	0313-6121-09	0303-6121-09	0312-6121-09	0302-6121-09	0314-6121-09	0304-6121-09
Orange	0317-6121-10	0307-6121-10	0313-6121-10	0303-6121-10	0312-6121-10	0302-6121-10	0314-6121-10	0304-6121-10
Pink	0317-6121-12	0307-6121-12	0313-6121-12	0303-6121-12	0312-6121-12	0302-6121-12	0314-6121-12	0304-6121-12
Tan	0317-6121-15	0307-6121-15	0313-6121-15	0303-6121-15	0312-6121-15	0302-6121-15	0314-6121-15	0304-6121-15
White	0317-6121-16	0307-6121-16	0313-6121-16	0303-6121-16	0312-6121-16	0302-6121-16	0314-6121-16	0304-6121-16
Yellow	0317-6121-17	0307-6121-17	0313-6121-17	0303-6121-17	0312-6121-17	0302-6121-17	0314-6121-17	0304-6121-17

PrintAID Microscope Slides



Dimensions: 25.0mm×75.0mm, Thickness: 1.0-1.2mm

Recommended for marking with the conventional labeling system and the thermal transfer slide labelers, including Thermo Slide-Mate® (Thermo Fisher) slide printer, and Signature® (Primer) slide printers

- Color frosted area of approx. 20 mm depth displays at one surface end
- Color frosted area is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Various colors are supplied



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-7101-02	0303-7101-02	0312-7101-02	0302-7101-02
● Green	0313-7101-08	0303-7101-08	0312-7101-08	0302-7101-08
● Purple	0313-7101-09	0303-7101-09	0312-7101-09	0302-7101-09
● Orange	0313-7101-10	0303-7101-10	0312-7101-10	0302-7101-10
● Pink	0313-7101-12	0303-7101-12	0312-7101-12	0302-7101-12
● Tan	0313-7101-15	0303-7101-15	0312-7101-15	0302-7101-15
○ White	0313-7101-16	0303-7101-16	0312-7101-16	0302-7101-16
● Yellow	0313-7101-17	0303-7101-17	0312-7101-17	0302-7101-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-7121-02	0303-7121-02	0312-7121-02	0302-7121-02
● Green	0313-7121-08	0303-7121-08	0312-7121-08	0302-7121-08
● Purple	0313-7121-09	0303-7121-09	0312-7121-09	0302-7121-09
● Orange	0313-7121-10	0303-7121-10	0312-7121-10	0302-7121-10
● Pink	0313-7121-12	0303-7121-12	0312-7121-12	0302-7121-12
● Tan	0313-7121-15	0303-7121-15	0312-7121-15	0302-7121-15
○ White	0313-7121-16	0303-7121-16	0312-7121-16	0302-7121-16
● Yellow	0313-7121-17	0303-7121-17	0312-7121-17	0302-7121-17


PrintAID Microscope Slides



Dimensions: 25.0mm×75.0mm, Thickness: 1.0-1.2mm

Recommended for marking with the conventional labeling system and the inkjet slide labelers, special for Inkjet IPS (Leica) slide printers

- Color frosted area of approx. 20 mm depth displays at one surface end
- Color frosted area is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Various colors are supplied



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-4101-02	0303-4101-02	0312-4101-02	0302-4101-02
● Green	0313-4101-08	0303-4101-08	0312-4101-08	0302-4101-08
● Purple	0313-4101-09	0303-4101-09	0312-4101-09	0302-4101-09
● Orange	0313-4101-10	0303-4101-10	0312-4101-10	0302-4101-10
● Pink	0313-4101-12	0303-4101-12	0312-4101-12	0302-4101-12
● Tan	0313-4101-15	0303-4101-15	0312-4101-15	0302-4101-15
○ White	0313-4101-16	0303-4101-16	0312-4101-16	0302-4101-16
● Yellow	0313-4101-17	0303-4101-17	0312-4101-17	0302-4101-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-4121-02	0303-4121-02	0312-4121-02	0302-4121-02
● Green	0313-4121-08	0303-4121-08	0312-4121-08	0302-4121-08
● Purple	0313-4121-09	0303-4121-09	0312-4121-09	0302-4121-09
● Orange	0313-4121-10	0303-4121-10	0312-4121-10	0302-4121-10
● Pink	0313-4121-12	0303-4121-12	0312-4121-12	0302-4121-12
● Tan	0313-4121-15	0303-4121-15	0312-4121-15	0302-4121-15
○ White	0313-4121-16	0303-4121-16	0312-4121-16	0302-4121-16
● Yellow	0313-4121-17	0303-4121-17	0312-4121-17	0302-4121-17


Photosensitive Microscope Slide



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

Recommended for marking with the conventional labeling system and the laser slide labelers.

- Color frosted area of approx. 20 mm depth displays at one surface end
- Color frosted area is resistant to all common chemicals and routine stains that are utilized in the laboratory
- Various colors are supplied



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-8101-02	0303-8101-02	0312-8101-02	0302-8101-02
● Green	0313-8101-08	0303-8101-08	0312-8101-08	0302-8101-08
● Purple	0313-8101-09	0303-8101-09	0312-8101-09	0302-8101-09
● Orange	0313-8101-10	0303-8101-10	0312-8101-10	0302-8101-10
● Pink	0313-8101-12	0303-8101-12	0312-8101-12	0302-8101-12
● Tan	0313-8101-15	0303-8101-15	0312-8101-15	0302-8101-15
○ White	0313-8101-16	0303-8101-16	0312-8101-16	0302-8101-16
● Yellow	0313-8101-17	0303-8101-17	0312-8101-17	0302-8101-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-8121-02	0303-8121-02	0312-8121-02	0302-8121-02
● Green	0313-8121-08	0303-8121-08	0312-8121-08	0302-8121-08
● Purple	0313-8121-09	0303-8121-09	0312-8121-09	0302-8121-09
● Orange	0313-8121-10	0303-8121-10	0312-8121-10	0302-8121-10
● Pink	0313-8121-12	0303-8121-12	0312-8121-12	0302-8121-12
● Tan	0313-8121-15	0303-8121-15	0312-8121-15	0302-8121-15
○ White	0313-8121-16	0303-8121-16	0312-8121-16	0302-8121-16
● Yellow	0313-8121-17	0303-8121-17	0312-8121-17	0302-8121-17


ColorInx™ Microscope Slides adaptable to Etching Slide Printers



Dimensions: 25.0mm×75.0mm, Thickness: 1.0-1.2mm

Recommended for marking with etching slide printers

- Color frosted areas of approx. 20 mm depth, both sides but at the same slide end for marking
- A black background color helps increase contrast on the marking area
- Color frosted areas are resistant to all common chemicals and routine stains that are utilized in the laboratory
- Various colors are supplied



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-5201-02	0303-5201-02	0312-5201-02	0302-5201-02
● Green	0313-5201-08	0303-5201-08	0312-5201-08	0302-5201-08
● Purple	0313-5201-09	0303-5201-09	0312-5201-09	0302-5201-09
● Orange	0313-5201-10	0303-5201-10	0312-5201-10	0302-5201-10
● Pink	0313-5201-12	0303-5201-12	0312-5201-12	0302-5201-12
● Tan	0313-5201-15	0303-5201-15	0312-5201-15	0302-5201-15
○ White	0313-5201-16	0303-5201-16	0312-5201-16	0302-5201-16
● Yellow	0313-5201-17	0303-5201-17	0312-5201-17	0302-5201-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-2221-02	0303-2221-02	0312-2221-02	0302-2221-02
● Green	0313-2221-08	0303-2221-08	0312-2221-08	0302-2221-08
● Purple	0313-2221-09	0303-2221-09	0312-2221-09	0302-2221-09
● Orange	0313-2221-10	0303-2221-10	0312-2221-10	0302-2221-10
● Pink	0313-2221-12	0303-2221-12	0312-2221-12	0302-2221-12
● Tan	0313-2221-15	0303-2221-15	0312-2221-15	0302-2221-15
○ White	0313-2221-16	0303-2221-16	0312-2221-16	0302-2221-16
● Yellow	0313-2221-17	0303-2221-17	0312-2221-17	0302-2221-17



ADHESION MICROSCOPE SLIDES


Adhesion Microscope Slides

The continuous liquid-base in stains as well as chemicals and the frequent washing throughout the immunohistochemical staining workflow, especially the immediate washing after the antigen-retrieval heating do cause great risk of tissue section loss or cell loss. The surface of the CITOTEST® Color-PLUS™ microscope slides is coated with a permanent especially-formulated adhesion which reliably binds the tissue sections to the glass, considerably minimizing the possibility of tissue-section or cell loss.

The CITOTEST® Color-PLUS™ microscope slides are also highly recommended in routine H&E staining to secure the problematic tissue sections like brain and spinal cord specimens.

CAUTIONS

1. Store the products under dry and room temperature conditions.
2. Use the products as soon as possible when the package is open.
3. The color frosted layer may fall off the glass where caustic conditions are applied.



Recommended ★
Strongly recommended ★ ★

Applications

	PCB	PCL	PCI
HISTOLOGY AND PATHOLOGY			
Paraffin embedded section	★ ★	★ ★	★ ★
Frozen section	★ ★	★ ★	★ ★
Liquid-rich specimens	★	★	★ ★
AUTOMATED IHC STAINER			
Roche Ventana BenchMark	★	★	★ ★
Leica BOND-MAX	★ ★	★ ★	★ ★
Dako Autostainer Link 48	★ ★	★ ★	★ ★
HIER PROCESSES			
Acid Retrieval Solution	★ ★	★ ★	★ ★
Alkaline Retrieval Solution	★	★	★ ★
CYTOLOGY			
LBC, membrane	★ ★	★ ★	★ ★
LBC, centrifugation	★ ★	★ ★	★ ★
Microbiology	★	★	★
Hematology	★	★	★

Color-PLUS™ PCB Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCB series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End which is resistant to most common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils and mark pens.

Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS		BEVELED EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7141-02	0303-7141-02	0312-7141-02	0302-7141-02	0314-7141-02	0304-7141-02
Green	0313-7141-08	0303-7141-08	0312-7141-08	0302-7141-08	0314-7141-08	0304-7141-08
Purple	0313-7141-09	0303-7141-09	0312-7141-09	0302-7141-09	0314-7141-09	0304-7141-09
Orange	0313-7141-10	0303-7141-10	0312-7141-10	0302-7141-10	0314-7141-10	0304-7141-10
Pink	0313-7141-12	0303-7141-12	0312-7141-12	0302-7141-12	0314-7141-12	0304-7141-12
Tan	0313-7141-15	0303-7141-15	0312-7141-15	0302-7141-15	0314-7141-15	0304-7141-15
White	0313-7141-16	0303-7141-16	0312-7141-16	0302-7141-16	0314-7141-16	0304-7141-16
Yellow	0313-7141-17	0303-7141-17	0312-7141-17	0302-7141-17	0314-7141-17	0304-7141-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS		BEVELED EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7151-02	0303-7151-02	0312-7151-02	0302-7151-02	0314-7151-02	0304-7151-02
Green	0313-7151-08	0303-7151-08	0312-7151-08	0302-7151-08	0314-7151-08	0304-7151-08
Purple	0313-7151-09	0303-7151-09	0312-7151-09	0302-7151-09	0314-7151-09	0304-7151-09
Orange	0313-7151-10	0303-7151-10	0312-7151-10	0302-7151-10	0314-7151-10	0304-7151-10
Pink	0313-7151-12	0303-7151-12	0312-7151-12	0302-7151-12	0314-7151-12	0304-7151-12
Tan	0313-7151-15	0303-7151-15	0312-7151-15	0302-7151-15	0314-7151-15	0304-7151-15
White	0313-7151-16	0303-7151-16	0312-7151-16	0302-7151-16	0314-7151-16	0304-7151-16
Yellow	0313-7151-17	0303-7151-17	0312-7151-17	0302-7151-17	0314-7151-17	0304-7151-17

Color-PLUS™ PCL Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCL series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils, mark pens as well as most thermal transfer based slide printer, including Thermo Slide-Mate® (Thermo Fisher) slide printer, and Signature® (Primer) slide printers

Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7181-02	0303-7181-02	0312-7181-02	0302-7181-02
Green	0313-7181-08	0303-7181-08	0312-7181-08	0302-7181-08
Purple	0313-7181-09	0303-7181-09	0312-7181-09	0302-7181-09
Orange	0313-7181-10	0303-7181-10	0312-7181-10	0302-7181-10
Pink	0313-7181-12	0303-7181-12	0312-7181-12	0302-7181-12
Tan	0313-7181-15	0303-7181-15	0312-7181-15	0302-7181-15
White	0313-7181-16	0303-7181-16	0312-7181-16	0302-7181-16
Yellow	0313-7181-17	0303-7181-17	0312-7181-17	0302-7181-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7191-02	0303-7191-02	0312-7191-02	0302-7191-02
Green	0313-7191-08	0303-7191-08	0312-7191-08	0302-7191-08
Purple	0313-7191-09	0303-7191-09	0312-7191-09	0302-7191-09
Orange	0313-7191-10	0303-7191-10	0312-7191-10	0302-7191-10
Pink	0313-7191-12	0303-7191-12	0312-7191-12	0302-7191-12
Tan	0313-7191-15	0303-7191-15	0312-7191-15	0302-7191-15
White	0313-7191-16	0303-7191-16	0312-7191-16	0302-7191-16
Yellow	0313-7191-17	0303-7191-17	0312-7191-17	0302-7191-17

Color-PLUS™ PCL Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCL series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils, mark pens as well as most Inkjet based slide printer, special for Inkjet IPS (Leica) slide printers
Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4181-02	0303-4181-02	0312-4181-02	0302-4181-02
Green	0313-4181-08	0303-4181-08	0312-4181-08	0302-4181-08
Purple	0313-4181-09	0303-4181-09	0312-4181-09	0302-4181-09
Orange	0313-4181-10	0303-4181-10	0312-4181-10	0302-4181-10
Pink	0313-4181-12	0303-4181-12	0312-4181-12	0302-4181-12
Tan	0313-4181-15	0303-4181-15	0312-4181-15	0302-4181-15
White	0313-4181-16	0303-4181-16	0312-4181-16	0302-4181-16
Yellow	0313-4181-17	0303-4181-17	0312-4181-17	0302-4181-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4191-02	0303-4191-02	0312-4191-02	0302-4191-02
Green	0313-4191-08	0303-4191-08	0312-4191-08	0302-4191-08
Purple	0313-4191-09	0303-4191-09	0312-4191-09	0302-4191-09
Orange	0313-4191-10	0303-4191-10	0312-4191-10	0302-4191-10
Pink	0313-4191-12	0303-4191-12	0312-4191-12	0302-4191-12
Tan	0313-4191-15	0303-4191-15	0312-4191-15	0302-4191-15
White	0313-4191-16	0303-4191-16	0312-4191-16	0302-4191-16
Yellow	0313-4191-17	0303-4191-17	0312-4191-17	0302-4191-17


Color-PLUS™ PCL Photosensitive Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCL series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils, mark pens as well as most laser based slide printer
Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8181-02	0303-8181-02	0312-8181-02	0302-8181-02
Green	0313-8181-08	0303-8181-08	0312-8181-08	0302-8181-08
Purple	0313-8181-09	0303-8181-09	0312-8181-09	0302-8181-09
Orange	0313-8181-10	0303-8181-10	0312-8181-10	0302-8181-10
Pink	0313-8181-12	0303-8181-12	0312-8181-12	0302-8181-12
Tan	0313-8181-15	0303-8181-15	0312-8181-15	0302-8181-15
White	0313-8181-16	0303-8181-16	0312-8181-16	0302-8181-16
Yellow	0313-8181-17	0303-8181-17	0312-8181-17	0302-8181-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8191-02	0303-8191-02	0312-8191-02	0302-8191-02
Green	0313-8191-08	0303-8191-08	0312-8191-08	0302-8191-08
Purple	0313-8191-09	0303-8191-09	0312-8191-09	0302-8191-09
Orange	0313-8191-10	0303-8191-10	0312-8191-10	0302-8191-10
Pink	0313-8191-12	0303-8191-12	0312-8191-12	0302-8191-12
Tan	0313-8191-15	0303-8191-15	0312-8191-15	0302-8191-15
White	0313-8191-16	0303-8191-16	0312-8191-16	0302-8191-16
Yellow	0313-8191-17	0303-8191-17	0312-8191-17	0302-8191-17


Color-PLUS™ PCL Microscope Slides adaptable to Etching Slide Printers



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The CITOTEST adhesion microscope slides adaptable to etching slide printers always come with a black color coat for contrasting on the opposite side of the marking color-coat side, markable with the conventional labeling system, pencils and mark pens as well as most etching slide printers. The color coating is resistant to most common chemicals and routine stains that are utilized in laboratory.

Recommended application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3281-02	0303-3281-02	0312-3281-02	0302-3281-02
Green	0313-3281-08	0303-3281-08	0312-3281-08	0302-3281-08
Purple	0313-3281-09	0303-3281-09	0312-3281-09	0302-3281-09
Orange	0313-3281-10	0303-3281-10	0312-3281-10	0302-3281-10
Pink	0313-3281-12	0303-3281-12	0312-3281-12	0302-3281-12
Tan	0313-3281-15	0303-3281-15	0312-3281-15	0302-3281-15
White	0313-3281-16	0303-3281-16	0312-3281-16	0302-3281-16
Yellow	0313-3281-17	0303-3281-17	0312-3281-17	0302-3281-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3291-02	0303-3291-02	0312-3291-02	0302-3291-02
Green	0313-3291-08	0303-3291-08	0312-3291-08	0302-3291-08
Purple	0313-3291-09	0303-3291-09	0312-3291-09	0302-3291-09
Orange	0313-3291-10	0303-3291-10	0312-3291-10	0302-3291-10
Pink	0313-3291-12	0303-3291-12	0312-3291-12	0302-3291-12
Tan	0313-3291-15	0303-3291-15	0312-3291-15	0302-3291-15
White	0313-3291-16	0303-3291-16	0312-3291-16	0302-3291-16
Yellow	0313-3291-17	0303-3291-17	0312-3291-17	0302-3291-17

Color-PLUS™ PCI Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface. End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, pencils, mark pens as well as most thermal transfer based slide printers.

Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (privileged use), manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7161-02	0303-7161-02	0312-7161-02	0302-7161-02
Green	0313-7161-08	0303-7161-08	0312-7161-08	0302-7161-08
Purple	0313-7161-09	0303-7161-09	0312-7161-09	0302-7161-09
Orange	0313-7161-10	0303-7161-10	0312-7161-10	0302-7161-10
Pink	0313-7161-12	0303-7161-12	0312-7161-12	0302-7161-12
Tan	0313-7161-15	0303-7161-15	0312-7161-15	0302-7161-15
White	0313-7161-16	0303-7161-16	0312-7161-16	0302-7161-16
Yellow	0313-7161-17	0303-7161-17	0312-7161-17	0302-7161-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7171-02	0303-7171-02	0312-7171-02	0302-7171-02
Green	0313-7171-08	0303-7171-08	0312-7171-08	0302-7171-08
Purple	0313-7171-09	0303-7171-09	0312-7171-09	0302-7171-09
Orange	0313-7171-10	0303-7171-10	0312-7171-10	0302-7171-10
Pink	0313-7171-12	0303-7171-12	0312-7171-12	0302-7171-12
Tan	0313-7171-15	0303-7171-15	0312-7171-15	0302-7171-15
White	0313-7171-16	0303-7171-16	0312-7171-16	0302-7171-16
Yellow	0313-7171-17	0303-7171-17	0312-7171-17	0302-7171-17

Color-PLUS™ PCI Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, pencils, mark pens as well as most inkjet based slide printers.

Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (privileged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4161-02	0303-4161-02	0312-4161-02	0302-4161-02
Green	0313-4161-08	0303-4161-08	0312-4161-08	0302-4161-08
Purple	0313-4161-09	0303-4161-09	0312-4161-09	0302-4161-09
Orange	0313-4161-10	0303-4161-10	0312-4161-10	0302-4161-10
Pink	0313-4161-12	0303-4161-12	0312-4161-12	0302-4161-12
Tan	0313-4161-15	0303-4161-15	0312-4161-15	0302-4161-15
White	0313-4161-16	0303-4161-16	0312-4161-16	0302-4161-16
Yellow	0313-4161-17	0303-4161-17	0312-4161-17	0302-4161-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4171-02	0303-4171-02	0312-4171-02	0302-4171-02
Green	0313-4171-08	0303-4171-08	0312-4171-08	0302-4171-08
Purple	0313-4171-09	0303-4171-09	0312-4171-09	0302-4171-09
Orange	0313-4171-10	0303-4171-10	0312-4171-10	0302-4171-10
Pink	0313-4171-12	0303-4171-12	0312-4171-12	0302-4171-12
Tan	0313-4171-15	0303-4171-15	0312-4171-15	0302-4171-15
White	0313-4171-16	0303-4171-16	0312-4171-16	0302-4171-16
Yellow	0313-4171-17	0303-4171-17	0312-4171-17	0302-4171-17

Color-PLUS™ PCI Photosensitive Microscope Slides




Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, pencils, mark pens as well as most laser based slide printers.

Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (privileged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8161-02	0303-8161-02	0312-8161-02	0302-8161-02
Green	0313-8161-08	0303-8161-08	0312-8161-08	0302-8161-08
Purple	0313-8161-09	0303-8161-09	0312-8161-09	0302-8161-09
Orange	0313-8161-10	0303-8161-10	0312-8161-10	0302-8161-10
Pink	0313-8161-12	0303-8161-12	0312-8161-12	0302-8161-12
Tan	0313-8161-15	0303-8161-15	0312-8161-15	0302-8161-15
White	0313-8161-16	0303-8161-16	0312-8161-16	0302-8161-16
Yellow	0313-8161-17	0303-8161-17	0312-8161-17	0302-8161-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8171-02	0303-8171-02	0312-8171-02	0302-8171-02
Green	0313-8171-08	0303-8171-08	0312-8171-08	0302-8171-08
Purple	0313-8171-09	0303-8171-09	0312-8171-09	0302-8171-09
Orange	0313-8171-10	0303-8171-10	0312-8171-10	0302-8171-10
Pink	0313-8171-12	0303-8171-12	0312-8171-12	0302-8171-12
Tan	0313-8171-15	0303-8171-15	0312-8171-15	0302-8171-15
White	0313-8171-16	0303-8171-16	0312-8171-16	0302-8171-16
Yellow	0313-8171-17	0303-8171-17	0312-8171-17	0302-8171-17

Color-PLUS™ PCI Microscope Slides adaptable to Etching Slide Printers




Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface. End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, pencils, mark pens as well as most etching based slide printers.

Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (privileged use), manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-3261-02	0303-3261-02	0312-3261-02	0302-3261-02
● Green	0313-3261-08	0303-3261-08	0312-3261-08	0302-3261-08
● Purple	0313-3261-09	0303-3261-09	0312-3261-09	0302-3261-09
● Orange	0313-3261-10	0303-3261-10	0312-3261-10	0302-3261-10
● Pink	0313-3261-12	0303-3261-12	0312-3261-12	0302-3261-12
● Tan	0313-3261-15	0303-3261-15	0312-3261-15	0302-3261-15
○ White	0313-3261-16	0303-3261-16	0312-3261-16	0302-3261-16
● Yellow	0313-3261-17	0303-3261-17	0312-3261-17	0302-3261-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90° CORNERS		GROUND EDGES, 45° CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
● Blue	0313-3271-02	0303-3271-02	0312-3271-02	0302-3271-02
● Green	0313-3271-08	0303-3271-08	0312-3271-08	0302-3271-08
● Purple	0313-3271-09	0303-3271-09	0312-3271-09	0302-3271-09
● Orange	0313-3271-10	0303-3271-10	0312-3271-10	0302-3271-10
● Pink	0313-3271-12	0303-3271-12	0312-3271-12	0302-3271-12
● Tan	0313-3271-15	0303-3271-15	0312-3271-15	0302-3271-15
○ White	0313-3271-16	0303-3271-16	0312-3271-16	0302-3271-16
● Yellow	0313-3271-17	0303-3271-17	0312-3271-17	0302-3271-17

Liquid-Base Cytology Microscope Slides


Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

Supplied for liquid-base cytology, e.g., TCT & LCT slide preparation.

- Ground Edges, 45°clipped corners
- Printed Circle or Bracket optional
- Coated and non-coated available
- Frosted area of approx. 20mm depth, displays at one surface end

Packing:

50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
72pcs/box, 1,440pcs/carton



Microscope Slides with Printed Circles

Compatible with single or double Cellfunnels for use in cyto centrifuges. Recommended for Thermo Scientific™ Shandon™ EZ Cytofunnel™ or other brands cytofunnels which have the same design.

MARKING AREA	DESCRIPTION	50PCS		72PCS	
		CSW	SUPER GRADE	CSW	SUPER GRADE
A	White Single printed circle Dia.: 13mm	0370-0003	0370-0103	0370-0203-16	0370-0303-16
	White Single printed circle Dia.: 13mm, coated	0370-0001	0370-0101	0370-0201-16	0370-0301-16
B	White Double printed circles Dia.: 13mm	0370-0004	0370-0104	0370-0204-16	0370-0304-16
	White Double printed circles Dia.: 13mm, coated	0370-0002	0370-0102	0370-0202-16	0370-0302-16
C	White Double printed circles Dia.: 13mm	0370-0020	0370-0120	0370-0220-16	0370-0320-16
	White Double printed circles Dia.: 13mm, coated	0370-0021	0370-0121	0370-0221-16	0370-0321-16

Microscope Slides with Printed Bracket


Specifically designed for use with the ThinPrep® and other brand processors.

MARKING AREA	DESCRIPTION	50PCS		72PCS	
		CSW	SUPER GRADE	CSW	SUPER GRADE
D	White Printed Bracket Dia.: 36mm, coated	0370-0013-16	0370-0113-16	0370-0213-16	0370-0313-16

MICROSCOPE SLIDES FOR BLOOD SMEAR

The microscope slides are treated with our up-to-date technology on the automated cleaning workflow, featuring debris-free surface which provides ideal wettability for blood smear preparation.

Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm




Auto-Smear™ Microscope Slides

Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm



Recommended for automatic blood smear, including Sysmex fully automated slide-maker Sp1000i, COULTER SlideMaker & SlideStainer Hematology System LH 755, UniCel® DxH SlideMaker & SlideStainer Coulter Cellular Analysis System and COULTER® Hematology Automation

Markable with the conventional labeling system, pencil, mark pens as well as the most popular thermal transfer based slide labelers, such as the Thermo Slide-Mate® (Thermo Fisher) and Signature® (Primer) slide printers



- Color frosted area of approx. 20 mm depth, displays at one surface end
- Color frosted end is resistant to most common chemicals and routine stains that are utilized in the laboratory
- Cleanliness and wettability of slide surfaces are strictly and uniformly controlled

MARKING AREA	CSW	SUPER GRADE	PACKING
● Pink	5860-0001	5860-0101	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
● Pink	5860-0002	5860-0102	72pcs/box, 1,440pcs/carton
○ White	5860-0011	5860-0111	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
○ White	5860-0012	5860-0112	72pcs/box, 1,440pcs/carton
● Green	5860-0045	5860-0145	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
● Green	5860-0046	5860-0146	72pcs/box, 1,440pcs/carton




PapSmear™ Microscope Slides

Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

Featuring with ideal wetability, suitable for manual blood smear preparations.

Instruction for Blood Smear



1. Place a drop on one end of slide, Hold the spatula at a 30° ~45°angle, toward the drop of blood
2. Wait until the blood spreads along the entire width of the spatula
3. Hold the spatula at the same angle, push it forward rapidly and smoothly



MARKING AREA	GROUND EDGES, 45° CORNERS		BEVELED EDGES, 45° CORNERS		PACKING
	CSW	SUPER GRADE	CSW	SUPER GRADE	
White	5860-0015	5860-0115	5860-0017	5860-0117	50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton 72pcs/box, 1,440pcs/carton
White	5860-0016	5860-0116	5860-0018	5860-0118	

Plastic Spatula For Blood Smear

Disposables for even, smooth spreading of blood onto a microscopic slide

- PS material , light weight to prevent blood cells from being destroyed
- Tapered ends, to help achieve uniform smears
- Usable on both sides, to meet requirements with different smear areas
- No sharp edges or points for safety, to avoid hands and fingers from cutting

P/N	DIM. (L×W, mm)	PACKING
5814-0001	75×26	50pcs/box
5814-0002	75×26	72pcs/box



**DIAGNOSTIC
MICROSCOPE SLIDES**


Diagnostic Microscope Slides



CITOTEST is a leading professional manufacturer of Microscope slides and Coverslips in China

CITOTEST® slides have been exported to various fields required microscope all over the world and widely used in different kinds of lab every day. In order to satisfy more high-end professional clients' requirements, a new production of Diagnostic microscope slides is being launched, including Immunofluorescence assay slides, Quality control slides, Counting chamber slides, Milk Smear Slides, and Latex glass slides. These slides can be used in laboratories of Histology, Virology, Serology, Microbiology, etc.


- All slides are printed on the high quality of CITOTEST® "Pathology grade" slides.
- Super white glass chosen features optically clear and lower autofluorescence.
- Standard microscope slides have the dimension of 75.0x25.0x1.0mm.
- All slides have 90° corners and ground edges.
- Plain and frosted end are provided depending on the printed design.
- The available colors are Black, White, Blue, Green, Orange, and Red.



Immunofluorescence Assay Slides

Immunofluorescence assay (IFA) slides are designed for use in immunofluorescence tests. IFA slides have wells created with an extreme hydrophobic Teflon coating. The proprietary coating has a contact angle greater than 150° with distilled water and approx. 20 microns of thickness, which can effectively prevent crosscontamination and fully increase well capacity. While the wells are hydrophilic to help antigens spread and enhance cell attachment.

Different kinds of well pattern can be provided. The number of wells present ranges from single to 30, the diameter of minimum well is 3mm. Available colors are White, Black, Blue, Green and Orange. IFA slides are also available with epoxy, adhesive coatings and frosted end slides on custom-design requirements.



The adhesive coating features a permanent positive charge on the surface of slide, providing the slide a firm electrostatic attraction to bind tissue sections, cells and cytology preparations, without the risk of tissue and cells loss during frequent washing when carrying out staining. Adhesive microscope slides are ideal for routine histology, immunohistochemical (IHC) staining and in-situ hybridization (ISH) tests.

Teflon Coated Microscope slides

An extremely Hydrophobic coating with Teflon ink has a contact angle of 120° with distilled water, even higher. These slides have a matte finish with a coated ink thickness of 20-30 um. Teflon coated slides are ideal for use in immunohistochemical (IHC) and immunofluorescence (IFA) tests with histological sections or cellular cultures.

- Super white glass.
- Dim. 75.0x25.0x1.0mm.
- 90° corner and ground edges.
- Plain and a frosted end depended on the printed design.
- Available colors: Black, White, Blue, Green, Orange, and Red.

PLAIN	10MM FROSTED END	20MM FROSTED END	DESCRIPTION	PACKING
0383-0000-01	0383-0100-01	0383-0200-01	1 well, 12mm	100pcs/box,2000pcs/carton
0383-0001-01	0383-0101-01	0383-0201-01	1 well, 15mm	100pcs/box,2000pcs/carton
0383-0002-01	0383-0102-01	0383-0202-01	1 well, 19mm	100pcs/box,2000pcs/carton
0383-0003-01	0383-0103-01	0383-0203-01	1 Oblong, 17*24mm	100pcs/box,2000pcs/carton
0383-0004-01	0383-0104-01	0383-0204-01	2 wells, 15mm	100pcs/box,2000pcs/carton
0383-0005-01	0383-0105-01	0383-0205-01	2 wells, 19mm	100pcs/box,2000pcs/carton
0383-0006-01	0383-0106-01	0383-0206-01	3 wells, 10mm	100pcs/box,2000pcs/carton
0383-0007-01	0383-0107-01	0383-0207-01	4 wells, 11mm	100pcs/box,2000pcs/carton
0383-0008-01	0383-0108-01	0383-0208-01	5 wells, 8mm	100pcs/box,2000pcs/carton
0383-0009-01	0383-0109-01	0383-0209-01	6 wells, 5mm	100pcs/box,2000pcs/carton
0383-0010-01	0383-0110-01	0383-0210-01	8 wells, 5mm	100pcs/box,2000pcs/carton
0383-0011-01	0383-0111-01	0383-0211-01	10 wells, 6mm	100pcs/box,2000pcs/carton
0383-0012-01	0383-0112-01	0383-0212-01	12 wells, 5mm	100pcs/box,2000pcs/carton
0383-0013-01	0383-0113-01	0383-0213-01	21 wells, 4mm	100pcs/box,2000pcs/carton
0383-0014-01	0383-0114-01	0383-0214-01	30 wells, 3mm	100pcs/box,2000pcs/carton

Note: *The printed color without specification will be automatically black.

Epoxy Coated Microscope slides

A hydrophobic coating with Epoxy ink has a contact angle approx. 90°-100° with distilled water. These slides have a gloss finish with a coated ink thickness of 10-20 um. Epoxy coated slides are used for diagnostic immunohistochemical (IHC) and immunofluorescence (IFA) tests with histological sections or cellular cultures.

- Super white glass.
- Dim. 75.0x25.0x1.0mm.
- 90° corner and ground edges.
- Plain and a frosted end depended on the printed design.
- Available colors: Black, White, Blue, Green, Orange, and Red.

PLAIN	10MM FROSTED END	20MM FROSTED END	DESCRIPTION	PACKING
0383-5000-01	0383-5100-01	0383-5200-01	1 well, 12mm	100pcs/box,2000pcs/carton
0383-5001-01	0383-5101-01	0383-5201-01	1 well, 15mm	100pcs/box,2000pcs/carton
0383-5002-01	0383-5102-01	0383-5202-01	1 well, 19mm	100pcs/box,2000pcs/carton
0383-5003-01	0383-5103-01	0383-5203-01	1 Oblong, 17*24mm	100pcs/box,2000pcs/carton
0383-5004-01	0383-5104-01	0383-5204-01	2 wells, 15mm	100pcs/box,2000pcs/carton
0383-5005-01	0383-5105-01	0383-5205-01	2 wells, 19mm	100pcs/box,2000pcs/carton
0383-5006-01	0383-5106-01	0383-5206-01	3 wells, 10mm	100pcs/box,2000pcs/carton
0383-5007-01	0383-5107-01	0383-5207-01	4 wells, 11mm	100pcs/box,2000pcs/carton
0383-5008-01	0383-5108-01	0383-5208-01	5 wells, 8mm	100pcs/box,2000pcs/carton
0383-5009-01	0383-5109-01	0383-5209-01	6 wells, 5mm	100pcs/box,2000pcs/carton
0383-5010-01	0383-5110-01	0383-5210-01	8 wells, 5mm	100pcs/box,2000pcs/carton
0383-5011-01	0383-5111-01	0383-5211-01	10 wells, 6mm	100pcs/box,2000pcs/carton
0383-5012-01	0383-5112-01	0383-5212-01	12 wells, 5mm	100pcs/box,2000pcs/carton
0383-5013-01	0383-5113-01	0383-5213-01	21 wells, 4mm	100pcs/box,2000pcs/carton
0383-5014-01	0383-5114-01	0383-5214-01	30 wells, 3mm	100pcs/box,2000pcs/carton

Note: *The printed color without specification will be automatically black.



Adhesion Microscope Slides

Adhesive microscope slides are ideal for routine histology, immuno-histochemical (IHC) staining and in-situ hybridization (ISH) tests.

- Super white glass.
- Dim. 75.0x25.0x1.0mm.
- 90° corner and ground edges.
- Plain and a frosted end depended on the printed design.
- Available colors: Black, White, Blue, Green, Orange, and Red.

Teflon Coated Adhesion Microscope Slides

An extremely Hydrophobic coating with Teflon ink has a contact angle of 120° with distilled water, even higher. These slides have a matte finish with a coated ink thickness of 20-30 um. Teflon coated slides are ideal for use in immunohistochemical (IHC) and immunofluorescence (IFA) tests with histological sections or cellular cultures.

PLAIN	10MM FROSTED END	20MM FROSTED END	DESCRIPTION	PACKING
0383-2000-01	0383-2100-01	0383-2200-01	1 well, 12mm	100pcs/box,2000pcs/carton
0383-2001-01	0383-2101-01	0383-2201-01	1 well, 15mm	100pcs/box,2000pcs/carton
0383-2002-01	0383-2102-01	0383-2202-01	1 well, 19mm	100pcs/box,2000pcs/carton
0383-2003-01	0383-2103-01	0383-2203-01	1 Oblong, 17*24mm	100pcs/box,2000pcs/carton
0383-2004-01	0383-2104-01	0383-2204-01	2 wells, 15mm	100pcs/box,2000pcs/carton
0383-2005-01	0383-2105-01	0383-2205-01	2 wells, 19mm	100pcs/box,2000pcs/carton
0383-2006-01	0383-2106-01	0383-2206-01	3 wells, 10mm	100pcs/box,2000pcs/carton
0383-2007-01	0383-2107-01	0383-2207-01	4 wells, 11mm	100pcs/box,2000pcs/carton
0383-2008-01	0383-2108-01	0383-2208-01	5 wells, 8mm	100pcs/box,2000pcs/carton
0383-2009-01	0383-2109-01	0383-2209-01	6 wells, 5mm	100pcs/box,2000pcs/carton
0383-2010-01	0383-2110-01	0383-2210-01	8 wells, 5mm	100pcs/box,2000pcs/carton
0383-2011-01	0383-2111-01	0383-2211-01	10 wells, 6mm	100pcs/box,2000pcs/carton
0383-2012-01	0383-2112-01	0383-2212-01	12 wells, 5mm	100pcs/box,2000pcs/carton
0383-2013-01	0383-2113-01	0383-2213-01	21 wells, 4mm	100pcs/box,2000pcs/carton
0383-2014-01	0383-2114-01	0383-2214-01	30 wells, 3mm	100pcs/box,2000pcs/carton

Note: *The printed color without specification will be automatically black.

Epoxy coated Adhesion Microscope Slides

A hydrophobic coating with Epoxy ink has a contact angle approx. 90°-100° with distilled water. These slides have a gloss finish with a coated ink thickness of 10-20 um. Epoxy coated slides are used for diagnostic immunohistochemical (IHC) and immunofluorescence (IFA) tests with histological sections or cellular cultures.

PLAIN	10MM FROSTED END	20MM FROSTED END	DESCRIPTION	PACKING
0383-7000-01	0383-7100-01	0383-7200-01	1 well, 12mm	100pcs/box,2000pcs/carton
0383-7001-01	0383-7101-01	0383-7201-01	1 well, 15mm	100pcs/box,2000pcs/carton
0383-7002-01	0383-7102-01	0383-7202-01	1 well, 19mm	100pcs/box,2000pcs/carton
0383-7003-01	0383-7103-01	0383-7203-01	1 Oblong, 17*24mm	100pcs/box,2000pcs/carton
0383-7004-01	0383-7104-01	0383-7204-01	2 wells, 15mm	100pcs/box,2000pcs/carton
0383-7005-01	0383-7105-01	0383-7205-01	2 wells, 19mm	100pcs/box,2000pcs/carton
0383-7006-01	0383-7106-01	0383-7206-01	3 wells, 10mm	100pcs/box,2000pcs/carton
0383-7007-01	0383-7107-01	0383-7207-01	4 wells, 11mm	100pcs/box,2000pcs/carton
0383-7008-01	0383-7108-01	0383-7208-01	5 wells, 8mm	100pcs/box,2000pcs/carton
0383-7009-01	0383-7109-01	0383-7209-01	6 wells, 5mm	100pcs/box,2000pcs/carton
0383-7010-01	0383-7110-01	0383-7210-01	8 wells, 5mm	100pcs/box,2000pcs/carton
0383-7011-01	0383-7111-01	0383-7211-01	10 wells, 6mm	100pcs/box,2000pcs/carton
0383-7012-01	0383-7112-01	0383-7212-01	12 wells, 5mm	100pcs/box,2000pcs/carton
0383-7013-01	0383-7113-01	0383-7213-01	21 wells, 4mm	100pcs/box,2000pcs/carton
0383-7014-01	0383-7114-01	0383-7214-01	30 wells, 3mm	100pcs/box,2000pcs/carton

Note: *The printed color without specification will be automatically black.

Quality Control Slides




Quality Control Slides are designed with specific control areas. Tested samples and controls can be placed on the same slides, in favor of rapid and accurate staining procedures. The controls may be anything samples used in quality control, education and proficiency programs, including Histology, Bacteriology, Mycobacteriology, Parasitology, and/or Mycology. They can be distinguished by the shape and the color of control areas. Utilizing quality control slides can be used to verify performance of stain reagents, staining methods, and the performance of personnel interpreting microscopic examinations.

Histology Control Slide

Histology control slides are ideal for staining procedures that require a control. The control section on the slide is easily identified with a bright red or other color box. Patient tissue and controls can be placed on the same slide, helping to make positive identification. Slides are treated with a special adhesive which can hold the tissue sections firmly on the slide. The control box on back are applied for use in automated instruments.

- Super white glass.
- Dim. 75.0x25.0x1.0mm.
- 90° corner and ground edges.
- White marking area and Red control box.
- Adhesion coated surface optional.



Packing: 100pcs/box, 2000pcs/carton

UNCOATED	COATED	DESCRIPTION
0384-0000	0384-2000	Histology Control Slides, Patient/Control, with 20mm (3/4") of frosted white end, red box on front.
0384-0001	0384-2001	Histology Control Slides, Patient/Control, with 20mm (3/4") of frosted white end, red box on back.
0384-0002	0384-2002	Histology Control Slides, Patient/Control, with 25mm (1") of frosted white end, red box on back.

Microbiology Control Slides


Microbiology control slides are ideal for staining procedures that require a control. The control section on the slide is easily identified with a color circle, straight line or other patterns. The sample and controls can be placed on the same slide, helping to make positive identification. All slides can be treated with a special adhesive which can hold the cells firmly on the slide.



- Super white glass.
- Dim. 75.0x25.0x1.0mm.
- 90° corner and ground edges.
- White marking area and control patterns depended on requirements.
- Adhesion coated surface optional.

Packing: 100pcs/box, 2000pcs/carton

UNCOATED	COATED	DESCRIPTION
0384-5000	0384-7000	Gram Slides, with 20mm (3/4") of frosted white end, 13mm diameter of blue wells with (+) and (-) by wells.
0384-5001	0384-7001	Gram Slides, with 20mm (3/4") of frosted white end, 13mm diameter of blue wells with (+) by wells.
0384-5002	0384-7002	Gram Slides, with 20mm (3/4") of frosted white end, 13mm diameter of blue wells with (-) by wells.
0384-5003	0384-7003	AFB Slides, with 20mm (3/4") of frosted white end, 13mm diameter of red wells with (+) and (-) by wells.
0384-5004	0384-7004	AFB Slides, with 20mm (3/4") of frosted white end, 13mm diameter of red wells with (+) by wells.
0384-5005	0384-7005	AFB Slides, with 20mm (3/4") of frosted white end, 13mm diameter of red wells with (-) by wells.
0384-5006	0384-7006	Nocardia Slides, with 20mm (3/4") of frosted white end, 13mm diameter of orange wells with (+) and (-) by wells.
0384-5007	0384-7007	Mycology Slides, with 20mm (3/4") of frosted white end, 13mm diameter of black wells with (+) and (-) by wells.
0384-5008	0384-7008	Parasitology Slides, with 20mm (3/4") of frosted white end, 23.5mm of (+) and (-) green areas.



Counting Chamber Slides



The counting chamber slide consists of a standard microscope slide and a coverslip, is specially designed for sperm counting and motility analysis, it can be used for the concentration assessment of microalgae, cellular elements or any other cell suspension.

All slides are disposable, ready-to-use, avoiding any contamination or influence due to the re-use effect. When a small volume of sample is deposited on the open edge of the chamber the fluid is distributed homogeneously by capillary action. The special coating on slide is nontoxic to cells and the thickness of coating are precisely controlled, helping to form a monolayer of cell in favor of the observation and counting.

Sperm Counting Chamber Slides

The slide has several counting chambers depending on requirements, each with the same depth. 10 microns and 20 microns depth chambers are available.

• 10 microns depth chambers


For human and stallion sperm analysis as allow having all the spermatozoa focus in one layer.

They can be of 1, 2, 4 or 6 chambers.

• 20 microns depth chambers

The best configuration for most of the animal species. In formats of 2, 4, 6 or 8 chambers per slide.

P/N	DESCRIPTION	PACKING
2230-0100	Counting chamber slides, 1 chamber, 10 microns depth	25pcs/box
2230-0101	Counting chamber slides, 2 chamber, 10 microns depth	25pcs/box
2230-0102	Counting chamber slides, 4 chamber, 10 microns depth	25pcs/box
2230-0103	Counting chamber slides, 6 chamber, 10 microns depth	25pcs/box
2230-0200	Counting chamber slides, 2 chamber, 20 microns depth	25pcs/box
2230-0201	Counting chamber slides, 4 chamber, 20 microns depth	25pcs/box
2230-0202	Counting chamber slides, 8 chamber, 20 microns depth	25pcs/box
2230-0203	Counting chamber slides, 6 chamber, 20 microns depth	25pcs/box




Milk Smear Slides




These slides are special designed for counting of somatic cells in the milk. They feature hydrophobic coating which can keep the milk in the circle without any spreading, and provide the marking area for labeling. The diameter of circle on slide is 11.28mm, permits rapid smearing of 0.01ml of milk over an approx. 1cm². For easy to count and judge, Slides with indicator marks are also supplied on requirements. These new slides have horizontal and vertical diameter lines per circle, which will enable you quickly count along the diameter and improve the accuracy and reliability of detection.

- Super white glass.
- Dim. 75.0x25.0x1.0mm.
- 90° corner and ground edges.
- Indicator marks depended on requirements.
- Five circles, the colors of Blue, Red, Green and Yellow are optional.



P/N	DESCRIPTION	PACKING
2230-5000-02	Milk Smear slides, 5 circles, Blue hydrophobic coating	100pcs/box, 2000pcs/carton
2230-6000-02	Milk Smear slides, 5 circles, Blue hydrophobic coating, with indicator marks	100pcs/box, 2000pcs/carton
2230-5000-13	Milk Smear slides, 5 circles, Red hydrophobic coating	100pcs/box, 2000pcs/carton
2230-6000-13	Milk Smear slides, 5 circles, Red hydrophobic coating, with indicator marks	100pcs/box, 2000pcs/carton
2230-5000-08	Milk Smear slides, 5 circles, Green hydrophobic coating	100pcs/box, 2000pcs/carton
2230-6000-08	Milk Smear slides, 5 circles, Green hydrophobic coating, with indicator marks	100pcs/box, 2000pcs/carton
2230-5000-17	Milk Smear slides, 5 circles, Yellow hydrophobic coating	100pcs/box, 2000pcs/carton
2230-6000-17	Milk Smear slides, 5 circles, Yellow hydrophobic coating, with indicator marks	100pcs/box, 2000pcs/carton




Capillary Gap Slides




Supplied with a special adhesive technology, which can be used for processing the chemically bond fresh or formalin-fixed frozen tissue sections. Ideal for special stains, immunocytochemicals, etc

Color Frosted area of approx. 30mm depth, support for well-known brand printers, including Thermo Slide-Mate (Thermo Fisher) slide printer, and Signature (Primer) Slide printers.

Alternative thickness of color frosted area at 100µm or 75µm, Extra gap width allows slide pair to function as a standard pair of capillary-gap slides




- Super white glass.
- Dim. 75.0x25.0x1.0mm.
- 90° corner and ground edges.




P/N	DESCRIPTION	PACKING
0385-0000	Capillary Gap Slides, thickness of color frosted area at 75µm	100pcs/box, 2000pcs/carton
0385-0001	Capillary Gap Slides, thickness of color frosted area at 100µm	100pcs/box, 2000pcs/carton

Latex Glass Slides

Reusable glass slides for Latex tests.



The glass slide is specially designed for Latex Agglutination and Coagulation tests. The hydrophilic, glossy surface ensures specimens and reagents to be mixed thoroughly and efficiently, meanwhile the hydrophobic rings provide a clear demarcation of the active working area, to avoid the flow of reagents and specimens out of testing area during experimental operation. And the complete printed rear side in a sharp contrast with the active working area, provides an efficient agglutination platform to assist in both work and judgment.



- 2mm thick glass plate with ground edges, for optimal safety
- Color rings printed on the front side of the plates, to prevent liquid from over flow
- The plates with a black rear side, for use with white Latex
- The plates with a white rear side, for use with colored Latex

Applications

-Pregnancy testing
-Syphilis testing
-Rubella testing
- Rheumatoid Factors detection
-C Reactive Protein detection
-Anti Streptolysin O detection
-Anti DNP detection
-Heterophile antibody detection
-and more


Packing: 50pcs/box, 400pcs/carton

P/N	REAR SIDE	FRONT SIDE
0382-0001-13	Black	Red oval rings
0382-0001-02	Black	Blue oval rings
0382-0001-08	Black	Green oval rings
0382-0001-16	Black	White oval rings
0382-0002-13	White	Red oval rings
0382-0002-02	White	Blue oval rings
0382-0002-01	White	Black oval rings

Mega-Sized Microscope Slides

Supplied for use with large pathological tissue specimens
Packing: 25pcs/box, 100pcs/box

Ground Edges, 90°Square Corners, Thickness: 1.0-1.2mm



Regular Microscope Slides

C5W	SUPER GRADE	MARKING AREA	DIM.
0350-0001	0350-0101	Plain	76.0×52.0mm
0350-0002	0350-0102	Frosted at 2 ends of 1side, End with 7mm depth	76.0×52.0mm
0351-0001	0351-0101	Plain	76.0×39.0mm
0351-0002	0351-0102	Frosted at 2 ends of 1side, End with 7mm depth	76.0×39.0mm

Adhesion Microscope Slides

C5W	SUPER GRADE	MARKING AREA	DIM.
0350-2001	0350-2101	Plain	76.0×52.0mm
0350-2002	0350-2102	Frosted at 2 ends of 1side , End with 7mm depth	76.0×52.0mm
0351-2001	0351-2101	Plain	76.0×39.0mm
0351-2002	0351-2102	Frosted at 2 ends of 1side , End with 7mm depth	76.0×39.0mm

Microscope marking pen

The marker pen is ideal for labeling microscope slides, which widely used in histology, pathology and cytology.

- Ultra fine tip allows precious hand marking on almost any surface including glass and plastic.
- Permanent and waterproof ink is resistant to many laboratory chemicals and solvents in processing and staining.



P/N	DESCRIPTION	PACKING
6000-0008	Black, 1mm nib	12pcs/box, 180pcs/carton



CITOTEST manufactures and supplies a variety of microscope cover glass. Provided with different raw materials, personalized craftsmanships, various quality control standards and tailored packaging, the CITOTEST® microscope cover glass family covers the following categories:

All of EstaGlas™, BoroGlas™, EcoGlas™ and Circular microscope cover glass are supplied in three thickness formats for meeting up with different operational preference and demands, #1.0 (0.13mm-0.16mm), #1.5 (0.16mm-0.19mm) and #2.0 (0.19mm-0.22mm).



MICROSCOPE COVER GLASS

EstaGlas™ Microscope Cover Glass

Made of the CITOTEST super white (CSW) glass sheets, feature ideal optical properties, stable molecular structure, uniform flatness and uniform sizes.
Recommended application: manual workflow in histology, cytology, urinalysis and microbiology



Packing by Quantity

SPECIFICATION (mm)	THICKNESS#1.0 0.13mm-0.16mm	THICKNESS#1.5 0.16mm-0.19mm	THICKNESS#2.0 0.19mm-0.22mm	PACKING
18×18	0340-0130	0341-0130	0342-0130	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
20×20	0340-0630	0341-0630	0342-0630	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×22	0340-1130	0341-1130	0342-1130	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×24	0340-1630	0341-1630	0342-1630	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×32	0340-2710	0341-2710	0342-2710	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×32	0340-2810	0341-2810	0342-2810	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×40	0340-3110	0341-3110	0342-3110	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×40	0340-3210	0341-3210	0342-3210	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×50	0340-3510	0341-3510	0342-3510	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×50	0340-3610	0341-3610	0342-3610	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×55	0340-5010	0341-5010	0342-5010	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×60	0340-3910	0341-3910	0342-3910	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×60	0340-4010	0341-4010	0342-4010	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton

Packing by Weight

SPECIFICATION (mm)	THICKNESS#1.0 0.13mm-0.16mm	THICKNESS#1.5 0.16mm-0.19mm	THICKNESS#2.0 0.19mm-0.22mm	PACKING
18×18	0340-0150	0341-0150	0342-0150	1oz./box, 10oz./promo-pk, 200oz./carton
20×20	0340-0650	0341-0650	0342-0650	1oz./box, 10oz./promo-pk, 200oz./carton
22×22	0340-1150	0341-1150	0342-1150	1oz./box, 10oz./promo-pk, 200oz./carton
24×24	0340-1650	0341-1650	0342-1650	1oz./box, 10oz./promo-pk, 200oz./carton
22×32	0340-2750	0341-2750	0342-2750	1oz./box, 10oz./promo-pk, 200oz./carton
24×32	0340-2850	0341-2850	0342-2850	1oz./box, 10oz./promo-pk, 200oz./carton
22×40	0340-3150	0341-3150	0342-3150	1oz./box, 10oz./promo-pk, 200oz./carton
24×40	0340-3250	0341-3250	0342-3250	1oz./box, 10oz./promo-pk, 200oz./carton
22×50	0340-3550	0341-3550	0342-3550	1oz./box, 10oz./promo-pk, 200oz./carton
24×50	0340-3650	0341-3650	0342-3650	1oz./box, 10oz./promo-pk, 200oz./carton
24×55	0340-5050	0341-5050	0342-5050	1oz./box, 10oz./promo-pk, 200oz./carton
22×60	0340-3950	0341-3950	0342-3950	1oz./box, 10oz./promo-pk, 200oz./carton
24×60	0340-4050	0341-4050	0342-4050	1oz./box, 10oz./promo-pk, 200oz./carton

*EstaGlas™ microscope cover glass in special dimensions are also available upon request

BoroGlas™ Microscope Cover Glass

Made of pure white, chemically-resistant borosilicate 3.3 glass sheets, feature excellent optical properties, stable molecular structure, uniform flatness, uniform sizes, ideal flexibility as well as considerably extended shelf-life.
Recommended application: manual and automatic workflow in histology, cytology, urinalysis and microbiology



Packing by Quantity

SPECIFICATION (mm)	THICKNESS#1.0 0.13mm-0.16mm	THICKNESS#1.5 0.16mm-0.19mm	THICKNESS#2.0 0.19mm-0.22mm	PACKING
18×18	0330-0130	0331-0130	0332-0130	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
20×20	0330-0630	0331-0630	0332-0630	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×22	0330-1130	0331-1130	0332-1130	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×24	0330-1630	0331-1630	0332-1630	200pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×32	0330-2710	0331-2710	0332-2710	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×32	0330-2810	0331-2810	0332-2810	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×40	0330-3110	0331-3110	0332-3110	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×40	0330-3210	0331-3210	0332-3210	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×50	0330-3510	0331-3510	0332-3510	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×50	0330-3610	0331-3610	0332-3610	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×55	0330-5010	0331-5010	0332-5010	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
22×60	0330-3910	0331-3910	0332-3910	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton
24×60	0330-4010	0331-4010	0332-4010	100pcs/box, 1,000pcs/promo-pk, 20,000pcs/carton

Packing by Weight


SPECIFICATION (mm)	THICKNESS#1.0 0.13mm-0.16mm	THICKNESS#1.5 0.16mm-0.19mm	THICKNESS#2.0 0.19mm-0.22mm	PACKING
18×18	0330-0150	0331-0150	0332-0150	1oz./box, 10oz./promo-pk, 200oz./carton
20×20	0330-0650	0331-0650	0332-0650	1oz./box, 10oz./promo-pk, 200oz./carton
22×22	0330-1150	0331-1150	0332-1150	1oz./box, 10oz./promo-pk, 200oz./carton
24×24	0330-1650	0331-1650	0332-1650	1oz./box, 10oz./promo-pk, 200oz./carton
22×32	0330-2750	0331-2750	0332-2750	1oz./box, 10oz./promo-pk, 200oz./carton
24×32	0330-2850	0331-2850	0332-2850	1oz./box, 10oz./promo-pk, 200oz./carton
22×40	0330-3150	0331-3150	0332-3150	1oz./box, 10oz./promo-pk, 200oz./carton
24×40	0330-3250	0331-3250	0332-3250	1oz./box, 10oz./promo-pk, 200oz./carton
22×50	0330-3550	0331-3550	0332-3550	1oz./box, 10oz./promo-pk, 200oz./carton
24×50	0330-3650	0331-3650	0332-3650	1oz./box, 10oz./promo-pk, 200oz./carton
24×55	0330-5050	0331-5050	0332-5050	1oz./box, 10oz./promo-pk, 200oz./carton
22×60	0330-3950	0331-3950	0332-3950	1oz./box, 10oz./promo-pk, 200oz./carton
24×60	0330-4050	0331-4050	0332-4050	1oz./box, 10oz./promo-pk, 200oz./carton

* BoroGlas™ microscope cover glass in special dimensions are also available upon request

EcoGlas™ Microscope Cover Glass

Recommended for economical manual workflow in routine experiments

- Made of green soda-lime glass sheets
- Uniform flatness and uniform sizes



SPECIFICATION (mm)	THICKNESS#1.0 0.13mm-0.16mm	THICKNESS#1.5 0.16mm-0.19mm	THICKNESS#2.0 0.19mm-0.22mm	PACKING
18×18	0340-0120	0341-0120	0342-0120	100pcs/box, 1,000pcs/vacuum-pk, 50,000pcs/carton
20×20	0340-0620	0341-0620	0342-0620	100pcs/box, 1,000pcs/vacuum-pk, 50,000pcs/carton
22×22	0340-1120	0341-1120	0342-1120	100pcs/box, 1,000pcs/vacuum-pk, 50,000pcs/carton
24×24	0340-1620	0341-1620	0342-1620	100pcs/box, 1,000pcs/vacuum-pk, 50,000pcs/carton
22×32	0340-2720	0341-2720	0342-2720	100pcs/box, 1,000pcs/vacuum-pk, 50,000pcs/carton
24×32	0340-2820	0341-2820	0342-2820	100pcs/box, 1,000pcs/vacuum-pk, 50,000pcs/carton

SPECIFICATION (mm)	THICKNESS#1.0 0.13mm-0.16mm	THICKNESS#1.5 0.16mm-0.19mm	THICKNESS#2.0 0.19mm-0.22mm	PACKING
22×40	0340-3120	0341-3120	0342-3120	100pcs/box, 1,000pcs/vacuum-pk, 40,000pcs/carton
24×40	0340-3220	0341-3220	0342-3220	100pcs/box, 1,000pcs/vacuum-pk, 40,000pcs/carton

SPECIFICATION (mm)	THICKNESS#1.0 0.13mm-0.16mm	THICKNESS#1.5 0.16mm-0.19mm	THICKNESS#2.0 0.19mm-0.22mm	PACKING
22×50	0340-3520	0341-3520	0342-3520	100pcs/box, 1,000pcs/vacuum-pk, 30,000pcs/carton
24×50	0340-3620	0341-3620	0342-3620	100pcs/box, 1,000pcs/vacuum-pk, 30,000pcs/carton
24×55	0340-5020	0341-5020	0342-5020	100pcs/box, 1,000pcs/vacuum-pk, 30,000pcs/carton
22×60	0340-3920	0341-3920	0342-3920	100pcs/box, 1,000pcs/vacuum-pk, 30,000pcs/carton
24×60	0340-4020	0341-4020	0342-4020	100pcs/box, 1,000pcs/vacuum-pk, 30,000pcs/carton

* EcoGlas™ microscope cover glass in special dimensions are also available upon request

Circular Microscope Cover Glass

Recommended for manual workflow in histology, cytology, urinalysis and microbiology


- Made of the CITOTEST super white (CSW) glass sheets
- Ideal optical properties, stable molecular structure
- Uniform flatness and uniform sizes

Packing: 100pcs/box,1,000pcs/vacuum-pk, 50,000pcs/carton

DIMENSIONS (MM)	THICKNESS NR.1.0 0.13-0.16MM	THICKNESS NR.1.5 0.16-0.19MM	THICKNESS NR.2.0 0.19-0.22MM
Φ 8	0344-0820	0345-0820	0346-0820
Φ 10	0344-1020	0345-1020	0346-1020
Φ 11	0344-1120	0345-1120	0346-1120
Φ 12	0344-1220	0345-1220	0346-1220
Φ 13	0344-1320	0345-1320	0346-1320
Φ 14	0344-1420	0345-1420	0346-1420
Φ 15	0344-1520	0345-1520	0346-1520
Φ 16	0344-1620	0345-1620	0346-1620
Φ 18	0344-1820	0345-1820	0346-1820
Φ 20	0344-2020	0345-2020	0346-2020
Φ 22	0344-2220	0345-2220	0346-2220
Φ 24	0344-2420	0345-2420	0346-2420
Φ 25	0344-2520	0345-2520	0346-2520



* Circular microscope cover glass in special dimensions are also available upon request



Hemacytometer Cover Glass

Recommended for blood cell counting used with hemacytometer

- Made of the CITOTEST super grade glass sheets
- Thickness 0.4mm
- Ideal optical properties, stable molecular structure
- Uniform flatness and uniform size

Packing: 10pcs/pack,10pcs/inner box, 1,000packs/carton

Polished edges

DIMENSIONS(MM)	P/N	CORNER STYLE
20×26	5850-0001	Corners 90°
22×22	5850-0004	Corners 45°
24×24	5850-0005	Corners 45°

Cut edges

DIMENSIONS(MM)	P/N	CORNER STYLE
20×26	5850-0002	Corners 90°
22×22	5850-0003	Corners 45°



Non-Sticking Microscope Cover Glass

NEW

-ASM surface treatment

Borosilicate glass 3.3 material, with special surface treatment, especially suitable for automatic sealing machine


- High cleanliness, free cleaning, the surface is specially treated, and it is not easy to stick under normal storage conditions;
- Appropriate flexibility, withstand a certain degree of bending without breaking;
- The surface characteristics make the sealant evenly distributed to prevent air bubbles from entering, and the specimens can be stored for a long time.

Borosilicate 3.3 Glass

P/N	SPECIFICATION	PACKING
0330-3281	24×40mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0330-3681	24×50mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0330-4081	24×60mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0331-3281	24×40mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0331-3681	24×50mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0331-4081	24×60mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton

Super White Glass

P/N	SPECIFICATION	PACKING
0340-3281	24×40mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0340-3681	24×50mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0340-4081	24×60mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0341-3281	24×40mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0341-3681	24×50mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
0341-4081	24×60mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton



PROCESSING EMBEDDING CASSETTES

EM 101 Biopsy Processing/Embedding Cassettes

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With square holes in dimensions of 1.6mm × 1.6mm
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the cassette labelers

Packed in bulk by dispenser box: 250pcs/disp.box, 2,000pcs/carton
Tapped: 40pcs/stack, 200pcs/inner box, 2,000pcs/carton

COLOR	PACKED IN BULK BY DISPENSER BOX	TAPPED
○ White	0107-1100-16	0107-1130-16
● Yellow	0107-1100-17	0107-1130-17
● Green	0107-1100-08	0107-1130-08
● Pink	0107-1100-12	0107-1130-12
● Blue	0107-1100-02	0107-1130-02
● Tan	0107-1100-15	0107-1130-15
● Purple	0107-1100-09	0107-1130-09
● Peach	0107-1100-11	0107-1130-11
● Orange	0107-1100-10	0107-1130-10
● Aqua	0107-1100-00	0107-1130-00



EM 102 Tissue Processing/Embedding Cassettes

One-piece and snap-latch design allows convenient one-hand operations, and helps relock the lid firmly during the processing

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With rectangle holes in dimensions of 5.0mm × 1.0mm
- With two large labeling areas on the cassette sides and one slanted front writing surface

Packed in bulk by dispenser box: 250pcs/disp.box, 2,000pcs/carton
Tapped: 40pcs/stack, 200pcs/inner box, 2,000pcs/carton

45° anterior writing surface, adaptable to most of the cassette labelers

COLOR	PACKED IN BULK BY DISPENSER BOX	TAPPED
○ White	0106-1100-16	0106-1130-16
● Yellow	0106-1100-17	0106-1130-17
● Green	0106-1100-08	0106-1130-08
● Pink	0106-1100-12	0106-1130-12
● Blue	0106-1100-02	0106-1130-02
● Tan	0106-1100-15	0106-1130-15
● Purple	0106-1100-09	0106-1130-09
● Peach	0106-1100-11	0106-1130-11
● Orange	0106-1100-10	0106-1130-10
● Aqua	0106-1100-00	0106-1130-00





EM 104 Biopsy Processing/Embedding Cassettes


- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With square holes in dimensions of 0.9mm×0.9mm
- With two large labeling areas on the cassette sides and an improved slanted front writing surface of 45°angle, specially for the thermal transfer printing system

Packed in bulk by dispenser box

COLOR	P/N	PACKING
○ White	0109-1100-16	250pcs/disp.box, 2,000pcs/carton
● Yellow	0109-1100-17	250pcs/disp.box, 2,000pcs/carton
● Green	0109-1100-08	250pcs/disp.box, 2,000pcs/carton
● Pink	0109-1100-12	250pcs/disp.box, 2,000pcs/carton
● Blue	0109-1100-02	250pcs/disp.box, 2,000pcs/carton
● Tan	0109-1100-15	250pcs/disp.box, 2,000pcs/carton
● Purple	0109-1100-09	250pcs/disp.box, 2,000pcs/carton
● Peach	0109-1100-11	250pcs/disp.box, 2,000pcs/carton
● Orange	0109-1100-10	250pcs/disp.box, 2,000pcs/carton
● Aqua	0109-1100-00	250pcs/disp.box, 2,000pcs/carton

Tapped

COLOR	P/N	PACKING
○ White	0109-1130-16	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Yellow	0109-1130-17	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Green	0109-1130-08	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Pink	0109-1130-12	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Blue	0109-1130-02	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Tan	0109-1130-15	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Purple	0109-1130-09	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Peach	0109-1130-11	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Orange	0109-1130-10	40pcs/stack, 200pcs/inner box, 2,000pcs/carton
● Aqua	0109-1130-00	40pcs/stack, 200pcs/inner box, 2,000pcs/carton



EM 105 Biopsy Processing/Embedding Cassettes

Special hinge design locks the cassette and lid together firmly, and realizes more frequent opening & closing of the lid without any danger of specimen loss

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With square holes in dimensions of 0.9mm×0.9mm
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the cassette labelers
- Detachable lid easy to remove and re-insert
- pre-loaded by sleeves or tapes, considerably optimizing the user convenience



Unassembled Packed by dispenser box: 500cassettes/disp.box+500lids/disp.box, 2,000sets/carton

Assembled Packed by dispenser box: 500sets/disp.box, 2,000sets/carton

Pre-loaded In sleeves: 75cassettes/sleeve+ 500lids/disp.box, 3,000sets/carton

Tapped: 40cassettes/stack,400cassettes/inner box+ 500lids/disp.box, 2,000pcs/carton

COLOR	UNASSEMBLED PACKED BY DISPENSER BOX	ASSEMBLED PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0103-1100-16	0103-1102-16	0103-1103-16	0103-1132-16
● Yellow	0103-1100-17	0103-1102-17	0103-1103-17	0103-1132-17
● Green	0103-1100-08	0103-1102-08	0103-1103-08	0103-1132-08
● Pink	0103-1100-12	0103-1102-12	0103-1103-12	0103-1132-12
● Blue	0103-1100-02	0103-1102-02	0103-1103-02	0103-1132-02
● Tan	0103-1100-15	0103-1102-15	0103-1103-15	0103-1132-15
● Purple	0103-1100-09	0103-1102-09	0103-1103-09	0103-1132-09
● Peach	0103-1100-11	0103-1102-11	0103-1103-11	0103-1132-11
● Orange	0103-1100-10	0103-1102-10	0103-1103-10	0103-1132-10
● Aqua	0103-1100-00	0103-1102-00	0103-1103-00	0103-1132-00



* Cassette&Lid Unassembled



* Cassette&Lid Assembled

Unassembled packed by dispenser box: 500cassettes/disp.box+500lids/disp.box, 2,000pcs/carton
 Assembled packed by dispenser box: 500pcs/disp.box, 2,000pcs/carton
 Pre-loaded in sleeves: 75cassettes/sleeve+500lids/disp.box, 3,000pcs/carton
 Tapped: 40cassettes/stack,400cassettes/innerbox+500lids/disp.box, 2,000pcs/carton

COLOR	UNASSEMBLED PACKED BY DISPENSER BOX	ASSEMBLED PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0100-1100-16	0100-1102-16	0100-1103-16	0100-1132-16
● Yellow	0100-1100-17	0100-1102-17	0100-1103-17	0100-1132-17
● Green	0100-1100-08	0100-1102-08	0100-1103-08	0100-1132-08
● Pink	0100-1100-12	0100-1102-12	0100-1103-12	0100-1132-12
● Blue	0100-1100-02	0100-1102-02	0100-1103-02	0100-1132-02
● Tan	0100-1100-15	0100-1102-15	0100-1103-15	0100-1132-15
● Purple	0100-1100-09	0100-1102-09	0100-1103-09	0100-1132-09
● Peach	0100-1100-11	0100-1102-11	0100-1103-11	0100-1132-11
● Orange	0100-1100-10	0100-1102-10	0100-1103-10	0100-1132-10
● Aqua	0100-1100-00	0100-1102-00	0100-1103-00	0100-1132-00

EM 107 Tissue Processing/Embedding Cassettes

Special designed with much larger slanted front writing surface, provide simple one-direction snapping and inserting to lock the lid into the cassette, allow more frequent opening & closing of the lid, and easy & safe removal and inserting of lid from and into the cassette without any danger of specimen loss


- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With rectangle holes in dimensions of 5.0mm x 1.0mm
- With two large labeling areas on the cassette sides and one much larger slanted front writing surface of 30°angle, not adaptable to the cassette labelers

Ⓐ Packed in bulk by dispenser box

COLOR	P/N	PACKING
○ White	0106-1106-16	250pcs/disp.box, 2,000pcs/carton
● Yellow	0106-1106-17	250pcs/disp.box, 2,000pcs/carton
● Green	0106-1106-08	250pcs/disp.box, 2,000pcs/carton
● Pink	0106-1106-12	250pcs/disp.box, 2,000pcs/carton
● Blue	0106-1106-02	250pcs/disp.box, 2,000pcs/carton
● Tan	0106-1106-15	250pcs/disp.box, 2,000pcs/carton
● Purple	0106-1106-09	250pcs/disp.box, 2,000pcs/carton
● Peach	0106-1106-11	250pcs/disp.box, 2,000pcs/carton
● Orange	0106-1106-10	250pcs/disp.box, 2,000pcs/carton
● Aqua	0106-1106-00	250pcs/disp.box, 2,000pcs/carton

Ⓑ Assembled packed by dispenser box

COLOR	P/N	PACKING
○ White	0100-1107-16	500pcs/disp.box, 2,000pcs/carton
● Yellow	0100-1107-17	500pcs/disp.box, 2,000pcs/carton
● Green	0100-1107-08	500pcs/disp.box, 2,000pcs/carton
● Pink	0100-1107-12	500pcs/disp.box, 2,000pcs/carton
● Blue	0100-1107-02	500pcs/disp.box, 2,000pcs/carton
● Tan	0100-1107-15	500pcs/disp.box, 2,000pcs/carton
● Purple	0100-1107-09	500pcs/disp.box, 2,000pcs/carton
● Peach	0100-1107-11	500pcs/disp.box, 2,000pcs/carton
● Orange	0100-1107-10	500pcs/disp.box, 2,000pcs/carton
● Aqua	0100-1107-00	500pcs/disp.box, 2,000pcs/carton




EM120 Tissue Processing/Embedding Cassettes

Two-piece design provides simple inserting and fixing to lock the lid into the cassette securely, allows more frequent opening & closing of the lid, and easy & safe removal and inserting of lid from and into the cassette without any danger of specimen loss.
 Unique design of back-mounted tab, allows unobtrusive anterior marking area, and alternative opening method can support simultaneous identification of samples and sample information reading .

- Made of P.O.M material, 100% resistant to chemical action of histological solvents
- With rectangle holes in dimensions of 5.0mm × 1.0mm
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the cassette printers
- Pre-loaded by sleeves or tapes, considerably optimizing the user convenience

Assembled, packed by dispenser box: 500sets/disp.box, 2,000sets/carton
 Pre-loaded in sleeves: 75sets/sleeve, 3,000sets/carton
 Tapped: 40sets/stack, 400sets/innerbox, 4,000sets/carton

COLOR	ASSEMBLED, PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0123-1101-16	0123-1120-16	0123-1133-16
● Yellow	0123-1101-17	0123-1120-17	0123-1133-17
● Green	0123-1101-08	0123-1120-08	0123-1133-08
● Pink	0123-1101-12	0123-1120-12	0123-1133-12
● Blue	0123-1101-02	0123-1120-02	0123-1133-02
● Tan	0123-1101-15	0123-1120-15	0123-1133-15
● Purple	0123-1101-09	0123-1120-09	0123-1133-09
● Peach	0123-1101-11	0123-1120-11	0123-1133-11
● Orange	0123-1101-10	0123-1120-10	0123-1133-10
● Aqua	0123-1101-00	0123-1120-00	0123-1133-00



EM 111& EM121 Biopsy Processing/Embedding Cassettes

Two-piece design provides simple inserting and fixing to lock the lid into the cassette securely, allows more frequent opening & closing of the lid, and easy & safe removal and inserting of lid from and into the cassette without any danger of specimen loss.

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With square holes in dimensions of 0.9mm × 0.9mm
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the cassette labelers
- Pre-loaded by sleeves or tapes, considerably optimizing the user convenience


Assembled Packed by dispenser box: 500pcs/disp.box, 2,000pcs/carton
 Pre-loaded In sleeves: 75pcs/sleeve, 3,000pcs/carton
 Tapped: 40pcs/stack, 400pcs/inner box, 4,000pcs/carton

EM 111

COLOR	ASSEMBLED PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0103-1101-16	0103-1120-16	0103-1133-16
● Yellow	0103-1101-17	0103-1120-17	0103-1133-17
● Green	0103-1101-08	0103-1120-08	0103-1133-08
● Pink	0103-1101-12	0103-1120-12	0103-1133-12
● Blue	0103-1101-02	0103-1120-02	0103-1133-02
● Tan	0103-1101-15	0103-1120-15	0103-1133-15
● Purple	0103-1101-09	0103-1120-09	0103-1133-09
● Peach	0103-1101-11	0103-1120-11	0103-1133-11
● Orange	0103-1101-10	0103-1120-10	0103-1133-10
● Aqua	0103-1101-00	0103-1120-00	0103-1133-00

EM 121

COLOR	ASSEMBLED PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0124-1101-16	0124-1120-16	0124-1133-16
● Yellow	0124-1101-17	0124-1120-17	0124-1133-17
● Green	0124-1101-08	0124-1120-08	0124-1133-08
● Pink	0124-1101-12	0124-1120-12	0124-1133-12
● Blue	0124-1101-02	0124-1120-02	0124-1133-02
● Tan	0124-1101-15	0124-1120-15	0124-1133-15
● Purple	0124-1101-09	0124-1120-09	0124-1133-09
● Peach	0124-1101-11	0124-1120-11	0124-1133-11
● Orange	0124-1101-10	0124-1120-10	0124-1133-10
● Aqua	0124-1101-00	0124-1120-00	0124-1133-00



EM110 & EM112 Tissue Processing/Embedding Cassettes

Two-piece design provides simple inserting and fixing to lock the lid into the cassette securely, allows more frequent opening & closing of the lid, and easy & safe removal and inserting of lid from and into the cassette without any danger of specimen loss.
Special compact design which is considerably space-saving, helping in loading and stacking of more cassettes into the automated cassette labeler

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With rectangle holes in dimensions of 5.0mm × 1.0 mm
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the cassette labelers
- Pre-loaded by sleeves or tapes, considerably optimizing the user convenience

Assembled, packed by dispenser box: 500pcs/disp.box, 2,000pcs/carton
Pre-loaded in sleeves: 75pcs/sleeve, 3,000pcs/carton
Tapped: 40pcs/stack, 400pcs/innerbox, 4,000pcs/carton

EM110 Tissue Processing/Embedding Cassettes

One compartment, allowing for a greatly improved fluid exchange

COLOR	ASSEMBLED, PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0100-1101-16	0100-1120-16	0100-1133-16
● Yellow	0100-1101-17	0100-1120-17	0100-1133-17
● Green	0100-1101-08	0100-1120-08	0100-1133-08
● Pink	0100-1101-12	0100-1120-12	0100-1133-12
● Blue	0100-1101-02	0100-1120-02	0100-1133-02
● Tan	0100-1101-15	0100-1120-15	0100-1133-15
● Purple	0100-1101-09	0100-1120-09	0100-1133-09
● Peach	0100-1101-11	0100-1120-11	0100-1133-11
● Orange	0100-1101-10	0100-1120-10	0100-1133-10
● Aqua	0100-1101-00	0100-1120-00	0100-1133-00

EM112 Tissue Processing/Embedding Cassettes

Four compartments, for processing/ embedding four tissue specimens at the same time

COLOR	ASSEMBLED, PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0100-1400-16	0100-1420-16	0100-1433-16
● Yellow	0100-1400-17	0100-1420-17	0100-1433-17
● Green	0100-1400-08	0100-1420-08	0100-1433-08
● Pink	0100-1400-12	0100-1420-12	0100-1433-12
● Blue	0100-1400-02	0100-1420-02	0100-1433-02
● Tan	0100-1400-15	0100-1420-15	0100-1433-15
● Purple	0100-1400-09	0100-1420-09	0100-1433-09
● Peach	0100-1400-11	0100-1420-11	0100-1433-11
● Orange	0100-1400-10	0100-1420-10	0100-1433-10
● Aqua	0100-1400-00	0100-1420-00	0100-1433-00

EM 113 Tissue Processing/Embedding Cassettes

One-piece & snap-latch design allows convenient one-hand operations, and helps relock the lid firmly during the processing

35° anterior writing surface, compatible with Leica and Sakura cassette labelers.

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With rectangle holes in dimensions of 5.0mm × 1.0mm
- With two large labeling areas on the cassette sides and one slanted front writing surface

Packed in bulk by dispenser box: 250pcs/disp.box, 2,000pcs/carton
Tapped: 40pcs/stack, 200pcs/inner box, 2,000pcs/carton

COLOR	PACKED IN BULK BY DISPENSER BOX	TAPPED
○ White	0106-1112-16	0106-1137-16
● Yellow	0106-1112-17	0106-1137-17
● Green	0106-1112-08	0106-1137-08
● Pink	0106-1112-12	0106-1137-12
● Blue	0106-1112-02	0106-1137-02
● Tan	0106-1112-15	0106-1137-15
● Purple	0106-1112-09	0106-1137-09
● Peach	0106-1112-11	0106-1137-11
● Orange	0106-1112-10	0106-1137-10
● Aqua	0106-1112-00	0106-1137-00



EM 114 Tissue Processing/Embedding Cassettes

Two-piece design provides simple inserting and fixing to lock the lid into the cassette securely, allows more frequent opening and closing of the lid, and easy and safe removal and inserting of lid from and into the cassette without any danger of specimen loss.

Special compact design which is considerably space-saving, helping in loading and stacking of more cassettes into the automated cassette labeler


- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With rectangle holes in dimensions of 5.0mm × 1.0mm
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the cassette labelers
- Pre-loaded by sleeves or tapes, considerably optimizing the user convenience

Assembled, packed by dispenser box: 500pcs/disp.box, 2,000pcs/carton
Pre-loaded in sleeves: 75pcs/sleeve, 3,000pcs/carton
Tapped: 40pcs/stack, 400pcs/innerbox, 4,000pcs/carton



* Upgraded lock ensures safety during processing
* special compact design is considerable space-saving

COLOR	ASSEMBLED, PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0100-1115-16	0100-1122-16	0100-1134-16
● Yellow	0100-1115-17	0100-1122-17	0100-1134-17
● Green	0100-1115-08	0100-1122-08	0100-1134-08
● Pink	0100-1115-12	0100-1122-12	0100-1134-12
● Blue	0100-1115-02	0100-1122-02	0100-1134-02
● Tan	0100-1115-15	0100-1122-15	0100-1134-15
● Purple	0100-1115-09	0100-1122-09	0100-1134-09
● Peach	0100-1115-11	0100-1122-11	0100-1134-11
● Orange	0100-1115-10	0100-1122-10	0100-1134-10
● Aqua	0100-1115-00	0100-1122-00	0100-1134-00



EM 117 Tissue Processing/Embedding Cassettes


One-piece & snap-latch design allows convenient one-hand operations, and helps relock the lid firmly during the processing

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With rectangle holes in dimensions of 5.0mm × 1.0 mm
- With two large labeling areas on the cassette sides and one slanted front writing surface
- Allowing for a greatly improved fluid exchange with extra anterior & posterior slots

Improved slanted front writing surface of 45°angle, specially for the thermal transfer cassette labelers

Packed in bulk by dispenser box: 250pcs/disp.box, 2,000pcs/carton
Tapped: 40pcs/stack, 200pcs/inner box, 2,000pcs/carton

COLOR	PACKED IN BULK BY DISPENSER BOX	TAPPED
○ White	0106-1117-16	0106-1138-16
● Yellow	0106-1117-17	0106-1138-17
● Green	0106-1117-08	0106-1138-08
● Pink	0106-1117-12	0106-1138-12
● Blue	0106-1117-02	0106-1138-02
● Tan	0106-1117-15	0106-1138-15
● Purple	0106-1117-09	0106-1138-09
● Peach	0106-1117-11	0106-1138-11
● Orange	0106-1117-10	0106-1138-10
● Aqua	0106-1117-00	0106-1138-00



EM 118 Tissue Processing/Embedding Cassettes

Two-piece design provides simple inserting and fixing to lock the lid into the cassette securely, allows more frequent opening and closing of the lid, and easy & safe removal and inserting of lid from and into the cassette without any danger of specimen loss.

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With rectangle holes in dimensions of 5.0mm × 1.0 mm
- With two large labeling areas on the cassette sides and one slanted front writing surface
- Allowing for a greatly improved fluid exchange with extra anterior & posterior slots

Improved slanted front writing surface of 45°angle, specially for the thermal transfer cassette labelers

Packed in bulk by dispenser box: 500pcs/disp.box, 2,000pcs/carton
Tapped(1): 40cassettes/stack, 400cassettes/inner box + 500lids/disp.box 2,000pcs/carton
Tapped(2): 40pcs/stack, 400pcs/inner box, 4,000pcs/carton

COLOR	PACKED IN BULK BY DISPENSER BOX	TAPPED(1)	TAPPED(2)
○ White	0122-1105-16	0122-1132-16	0122-1133-16
● Yellow	0122-1105-17	0122-1132-17	0122-1133-17
● Green	0122-1105-08	0122-1132-08	0122-1133-08
● Pink	0122-1105-12	0122-1132-12	0122-1133-12
● Blue	0122-1105-02	0122-1132-02	0122-1133-02
● Tan	0122-1105-15	0122-1132-15	0122-1133-15
● Purple	0122-1105-09	0122-1132-09	0122-1133-09
● Peach	0122-1105-11	0122-1132-11	0122-1133-11
● Orange	0122-1105-10	0122-1132-10	0122-1133-10
● Aqua	0122-1105-00	0122-1132-00	0122-1133-00

EM 115&116 Biopsy Processing Cassettes

Two-piece design and features with large anterior & posterior slots, especially used for processing/embedding tiny biopsy specimens without using biopsy pad, locking the lid securely without any danger of specimen loss and achieving ideal experimental effect

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With 1,280 square holes in dimensions of 0.26mm × 0.26 mm, and over 40 large rectangle holes to maximize fluid exchange, thereby ensuring proper drainage
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the printers labelers
- pre-loaded by sleeves or tapes, considerably optimizing the user convenience

Assembled Packed by dispenser box: 500pcs/disp.box, 2,000pcs/carton
Pre-loaded In sleeves: 75pcs/sleeve, 3,000pcs/carton
Tapped: 40pcs/stack, 400pcs/inner box, 4,000pcs/carton

EM 115

One compartment, allowing for a greatly improved fluid exchange instead of using biopsy pad.

COLOR	ASSEMBLED, PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0121-1101-16	0121-1130-16	0121-1133-16
● Yellow	0121-1101-17	0121-1130-17	0121-1133-17
● Green	0121-1101-08	0121-1130-08	0121-1133-08
● Pink	0121-1101-12	0121-1130-12	0121-1133-12
● Blue	0121-1101-02	0121-1130-02	0121-1133-02
● Tan	0121-1101-15	0121-1130-15	0121-1133-15
● Purple	0121-1101-09	0121-1130-09	0121-1133-09
● Peach	0121-1101-11	0121-1130-11	0121-1133-11
● Orange	0121-1101-10	0121-1130-10	0121-1133-10
● Aqua	0121-1101-00	0121-1130-00	0121-1133-00

EM 116

Four compartments, for processing/embedding four tiny biopsy specimens at the same time



COLOR	ASSEMBLED, PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0121-1401-16	0121-1430-16	0121-1433-16
● Yellow	0121-1401-17	0121-1430-17	0121-1433-17
● Green	0121-1401-08	0121-1430-08	0121-1433-08
● Pink	0121-1401-12	0121-1430-12	0121-1433-12
● Blue	0121-1401-02	0121-1430-02	0121-1433-02
● Tan	0121-1401-15	0121-1430-15	0121-1433-15
● Purple	0121-1401-09	0121-1430-09	0121-1433-09
● Peach	0121-1401-11	0121-1430-11	0121-1433-11
● Orange	0121-1401-10	0121-1430-10	0121-1433-10
● Aqua	0121-1401-00	0121-1430-00	0121-1433-00



Standard Tissue Processing/Embedding Cassettes

These models developed at the primary stage for specimen processing, embedding and sectioning, match with the CITOTEST metal lid#

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45°angle, adaptable to most of the cassette labelers
- Compact design for fitting securely in microtome chuck adapters, and preventing cassette deformation during the process of pathological section



Larger locking slot is designed to accept the metal lid.

EM 103 Tissue processing/Embedding cassettes

Made of P.O.M material, With round holes of Ø 2.4 mm

Bulk-packed by dispenser box: 500sets/disp.box, 2,000sets/carton

Pre-loaded in sleeves: 75pcs/sleeve, 3,000pcs/carton

Tapped: 40pcs/stack, 400pcs/inner box, 4,000pcs/carton

COLOR	BULK-PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0114-1104-16	0114-1123-16	0114-1135-16
● Yellow	0114-1104-17	0114-1123-17	0114-1135-17
● Green	0114-1104-08	0114-1123-08	0114-1135-08
● Pink	0114-1104-12	0114-1123-12	0114-1135-12
● Blue	0114-1104-02	0114-1123-02	0114-1135-02
● Tan	0114-1104-15	0114-1123-15	0114-1135-15
● Purple	0114-1104-09	0114-1123-09	0114-1135-09
● Peach	0114-1104-11	0114-1123-11	0114-1135-11
● Orange	0114-1104-10	0114-1123-10	0114-1135-10
● Aqua	0114-1104-00	0114-1123-00	0114-1135-00

EMSL 103 Stainless Steel Lids

Designed as a re-usable lid for the CITOTEST standard cassettes EM103 and EM108 and EM109, as well as most of the other standard cassettes in the market
Made of stainless steel material, with fine polished surface.

Packaging: 10pcs/pack, 250pcs/inner box, 2500pcs/carton

SPECIFICATION	P/N
Stainless steel cassette lid	0110-0005



EM108 Tissue Processing/Embedding Cassettes



Made of P.O.M material,
With rectangle holes in dimensions of 5.0mm×1.0mm

Bulk-packed by dispenser box: 500pcs/disp.box, 2,000pcs/carton

Pre-loaded in sleeves: 75pcs/sleeve, 3,000pcs/carton

Tapped: 40pcs/stack, 400pcs/inner box, 4,000pcs/carton

COLOR	BULK-PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0112-1100-16	0112-1123-16	0112-1135-16
● Yellow	0112-1100-17	0112-1123-17	0112-1135-17
● Green	0112-1100-08	0112-1123-08	0112-1135-08
● Pink	0112-1100-12	0112-1123-12	0112-1135-12
● Blue	0112-1100-02	0112-1123-02	0112-1135-02
● Tan	0112-1100-15	0112-1123-15	0112-1135-15
● Purple	0112-1100-09	0112-1123-09	0112-1135-09
● Peach	0112-1100-11	0112-1123-11	0112-1135-11
● Orange	0112-1100-10	0112-1123-10	0112-1135-10
● Aqua	0112-1100-00	0112-1123-00	0112-1135-00



EM 109 Tissue Processing/Embedding Cassettes


Made of P.O.M material,
With square holes in dimensions of 1.6mm×1.6mm

Bulk-packed by dispenser box: 500pcs/disp.box, 2,000pcs/carton

Pre-loaded in sleeves: 75pcs/sleeve, 3,000pcs/carton

Tapped: 40pcs/stack, 400pcs/inner box, 4,000pcs/carton

COLOR	BULK-PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
○ White	0113-1100-16	0113-1123-16	0113-1135-16
● Yellow	0113-1100-17	0113-1123-17	0113-1135-17
● Green	0113-1100-08	0113-1123-08	0113-1135-08
● Pink	0113-1100-12	0113-1123-12	0113-1135-12
● Blue	0113-1100-02	0113-1123-02	0113-1135-02
● Tan	0113-1100-15	0113-1123-15	0113-1135-15
● Purple	0113-1100-09	0113-1123-09	0113-1135-09
● Peach	0113-1100-11	0113-1123-11	0113-1135-11
● Orange	0113-1100-10	0113-1123-10	0113-1135-10
● Aqua	0113-1100-00	0113-1123-00	0113-1135-00



NEW

CITOTEST has newly launched a series of special processing/embedding cassettes, provided with 45° slanted angle which adapt to most laser cassette labeler in the market. Comparing to the inkjet labelers and thermal transfer labelers, the laser labeler does not require ink or impact ribbons, and the marking on cassettes can be permanent, indelible, especially meeting the most critical demands from laboratories from grossing to archiving, meanwhile, the marked cassettes can be scanned anywhere in the lab for quick and correct identification of block, in favor of improving workflow efficiencies and saving more time on diagnostics for lab technicians.

Processing/ Embedding Cassettes




Adapt to most of the laser cassettes labeler, providing reliable, permanent, indelible marking for specimen identification



Labeling parameters
 Laser: Optical fiber
 Wavelength: 1064nm
 Average power: 16W
 Hertz: 30kHz
 Speed: 1200m/sec

Made of P.O.M material, 100% resistant to chemical reaction of histological solvents; With one slanted front writing surface of 45° angle, adaptable to most of the laser cassette labelers; Pre-loaded by tapes or tubing dispensers, easy to be installed in the printing instruments; Available in ten colors



Pre-loaded In Sleeves

- (A) 75pcs/sleeve, 3,000pcs/carton
- (B) 75cassettes/sleeve + 500lids/disp.box, 3,000sets/carton

	P/N	COLOR
(A)	EM103 0114-2123-16	White
(B)	EM106 0100-2103-16	White
(B)	EM105 0103-2103-16	White
(A)	EM120 0123-2120-16	White
(A)	EM111 0103-2120-16	White

Tapped

- (A) 40pcs/stack, 400pcs/inner box, 4,000pcs/carton
- (B) 40pcs/stack, 200pcs/inner box, 2,000pcs/carton
- (C) 40pcs/stack, 400pcs/inner box +500lids/disp.box, 2,000sets/carton

	P/N	COLOR
(A)	EM103 0114-2135-16	White
(C)	EM106 0100-2132-16	White
(B)	EM102 0106-2130-16	White
(C)	EM105 0103-2132-16	White
(A)	EM120 0123-2133-16	White
(A)	EM111 0103-2133-16	White


Embedding Rings (EMOR201)

Designed for holding the tissue sample in the microtome chuck adapter for further sectioning

Packing: 500pcs/dispenser box, 2,000pcs/carton

ABS Material

COLOR	P/N
White	0260-0000-16
Yellow	0260-0000-17
Green	0260-0000-08
Pink	0260-0000-12
Blue	0260-0000-02
Tan	0260-0000-15
Purple	0260-0000-09
Peach	0260-0000-11
Orange	0260-0000-10
Aqua	0260-0000-00



Base Molds

Reusable Base Molds

Designed for facilitating same styles and same types of cassettes and embedding O rings in specimen embedding, compatible with most of the commonly-used embedding equipments

With fine polished surface and well-grounded corners, they are strong enough for limited re-use.



Packing: 10pcs/pack, 200pcs/inner box, 2,000pcs/carton

NET WELL DIM. (L×W)	EXTERNAL DIM. (L×W×H)	WITH NET WELL DEPTH OF 7MM	WITH NET WELL DEPTH OF 9MM FOR PROCESSING BIGGER TISSUE SPECIMENS
7×7	52×35×11	0203-0001	0203-1001
15×15	52×35×11	0203-0002	0203-1002
24×24	52×35×11	0203-0003	0203-1003
37×24	52×35×11	0203-0004	0203-1004
30×24	52×35×11	0203-0005	0203-1005

Disposable Base Molds

Designed for facilitating same styles and same types of cassettes and embedding "O" rings in specimen embedding, compatible with most of the commonly-used embedding equipments

With smooth surface finishing and round corners.
Specially-formulated plastic material



P/N	WELL DIM. (L×W)	EXTERNAL DIM. (L×W×H)	PACKING
0203-0006	7×7×5	50×37×12	500pcs/pack,5,000pcs/carton
0203-0007	17×17×5	50×37×12	500pcs/pack,5,000pcs/carton
0203-0008	25×25×5	50×37×12	500pcs/pack,5,000pcs/carton
0203-0009	25×32×5	50×37×12	500pcs/pack,5,000pcs/carton
0203-0010	25×40×5	50×37×12	500pcs/pack,5,000pcs/carton
0203-0011	30×20×12	55×38×18	50pcs/pack,500pcs/carton




Biopsy Bags

Used to hold small biopsy and histological specimens during processing

- Made of solvent resistant Polyester, effectively prevent the deformation
- Thin mesh (0.2×0.2 mm) provides excellent fluid exchange and avoids tiny specimen floss
- Convenient peel-open bag allows for placing and removing specimen rapidly

P/N	DIMENSION	PACKING
0110-0012	30×45 mm	200 pcs/box, 20,000 pcs/carton
0110-0013	45×75 mm	100 pcs/box, 20,000 pcs/carton
0110-0014	75×95 mm	100 pcs/box, 5,000 pcs/carton




1. Put the biopsy in a bag and fold its end
2. Place the bag into a cassette and close to process as a usual processing
3. Extract the bag from the cassette
4. Open it with care by extending the edges
5. extract the sample for embedding in paraffin

Biopsy Foam Pads

Designed for assisting the regular cassettes in holding the tiny biopsy specimens in place during processing, fits all the CITOTEST regular cassettes
Made of PU foam

Packing: 5,000pcs/pack, 20,000pcs/carton

P/N	DIMENSION (L×W×H)
0110-0001	30mm×25mm×2mm



Cassette Dispensing Sleeve Set

Used for pre-loading the cassettes, which can be accommodated directly by the cassettes labelers so as to eliminate hand inserting of cassettes and optimize the user convenience.

- With two sleeve end stoppers
- Accommodate up to approx.75pcs of the CITOTEST regular cassettes with lids

PVC Material

P/N	SPECIFICATION	PACKING
0110-0002	Cassette Dispensing Sleeve Set	200sleeves/carton+400 end stoppers/carton
0110-0003	Cassette Dispensing Sleeve	200sleeves/carton
0110-0004	End Stopper	850pcs/box,8x850pcs/carton

Paraffin Wax

Designed for improved tissue section infiltration and embedding and supplied in convenient zip-lock bag

- Compound of purified paraffin and plastic polymers
- Without dimethyl sulfoxide (DMSO)
- Delivers smooth cutting, minimized compression & cell disruption
- Cutting thickness lowered to 2 μm
- Convenient pearl shape
- Wide supply range from 52°C to 62°C for your choice



Packing: 1kg/bag, 10kg/carton

MELTING POINT °C	56-58	58-60	60-62	52-54	54-56
P/N	0200-0007	0200-0008	0200-0009	0200-0013	0200-0014

Microtome Blades

Supplied as disposables for facilitating the microtome and cryostat in tissue specimen sectioning, using the laterally displaceable blade holder, the quality microtome blades we supply below can be used as several blades just by the convenient cutting point re-positioning.

Stainless steel material



Made in Japan



Made in Germany

Packaging: 50pcs/box, 500pcs/inner box, 5000pcs/carton

P/N	BRAND	SUPPLIER PARTNO.
0390-0001	Feather	R35
0390-0002	Feather	A35
0390-0003	Feather	S35
0390-0004	Feather	N35
0391-0001	Leica	819
0391-0002	Leica	818



Staining Rack

Suitable for placing and storing microscope slides, accept both the CITOTEST standard microscope slides in dimensions of 75.0mm×25.0mm×(0.8-1.2)mm and the microscope slides of 76.0mm×26.0mm×(0.8-1.2)mm in the market

Stainless steel material

P/N	SPECIFICATION	PACKING
Ⓐ 0420-0001	20-PC, Used with the CITOTEST slide staining jar #0410-0004	1pcs/box, 50pcs/carton
Ⓑ 0420-0002	30-PC, Used with the CITOTEST slide staining jar #0410-0003	1pcs/box, 50pcs/carton
Ⓒ 0420-0003	60-PC, Used with the CITOTEST slide staining jar #0410-0005	1pcs/box, 50pcs/carton

Staining Jar

Suitable for microscope slides staining in routine operations, accept both the CITOTEST standard microscope slides in dimensions of 75.0mm×25.0mm×(0.8-1.2)mm and the microscope slides of 76.0mm×26.0mm×(0.8-1.2)mm on the market.



Mono™ Slide Staining Jar

- Made of borosilicate glass material
- With built-in slots for slide staining free from the use of separate slide racks

P/N	SPECIFICATION	PACKING
Ⓐ 0410-0000	5-PC, with glass lid, coplin type	1pcs/box, 20pcs/carton
Ⓑ 0410-0001	9-PC, with glass lid, hellendhal type	1pcs/box, 20pcs/carton
Ⓒ 0410-0002	10-PC, with glass lid, schifferdecker type	1pcs/box, 50pcs/carton



Sono™ Slide Staining Jar

- Made of borosilicate glass material
- Designed for slide staining together with related


P/N	SPECIFICATION	PACKING
Ⓐ 0410-0003	30-PC, with glass lid. Used with the CITOTEST slide staining rack #0420-0002	1pcs/box, 20pcs/carton
Ⓑ 0410-0004	20-PC, with glass lid, Papanicolau type Used with the CITOTEST slide staining rack #0420-0001	1pcs/box, 20pcs/carton
Ⓒ 0410-0005	60-PC, with glass lid. Used with the CITOTEST slide staining rack #0420-0003	1pcs/box, 10pcs/carton




Cardboard Tray, Mailer & Wallets

Designed as a convenient solution for temporary holding and transferring of microscope slides
Accommodate both the CITOTEST standard microscope slides in dimensions of 75.0mm×25.0mm×(0.8-1.2mm) and the microscope slides of 76.0mm×26.0mm×(0.8-1.2)mm in the market.

Made of heavy duty cardboard, reusable



P/N	SPECIFICATION	PACKING
(A) 0500-8001	1-pc, with lid and dividers	50pcs/inner box, 1,000pcs/carton
(B) 0500-8002	2-pc, with lid and dividers	50pcs/inner box, 500pcs/carton
(C) 0500-8003	3-pc, with lid and dividers	50pcs/inner box, 500pcs/carton
(D) 0500-8004	4-pc, with lid and dividers	50pcs/inner box, 200pcs/carton
(E) 0500-7009	9-pc, without lid or divider	10pcs/inner box, 250pcs/carton
(F) 0500-8010	10-pc, with lid and dividers	5pcs/inner box, 50pcs/carton
(G) 0500-8020	20-pc, with lid and dividers	5pcs/inner box, 50pcs/carton
(H) 0500-7020	20-pc, without lid, with dividers	10pcs/inner box, 100pcs/carton



P/N	SPECIFICATION	PACKING
(I) 0500-8030	30-pc, with lid and dividers	5pcs/inner box, 50pcs/carton
(J) 0500-8048	48-pc, with lid and dividers	5pcs/inner box, 50pcs/carton
(K) 0500-8024	24-pc, with lid, without divider	5pcs/inner box, 50pcs/carton
(L) 0500-8032	32-pc, with lid, without divider	5pcs/inner box, 50pcs/carton

**TRANSPORT
STORAGE
ARCHIVING**

Plastic Tray, Mailer & Wallets

Designed as a convenient solution for temporary holding and transferring of microscope slides
Accommodate both the CITOTEST standard microscope slides in dimensions of 75.0mm×25.0mm×(0.8-1.2)mm and the microscope slides of 76.0mm×26.0mm×(0.8-1.2)mm in the market.


Made of plastic material , reusable




P/N	SPECIFICATION	PACKING
Ⓐ 0500-1001	1-pc, NEUTRAL color. PP material	200pcs/dispenser box, 2,000pcs/carton
Ⓑ 0500-2002	2-pc, NEUTRAL color. PP material	100pcs/dispenser box, 1,000pcs/carton
Ⓒ 0500-1003	3-pc, NEUTRAL color. PP material	100pcs/dispenser box, 1,000pcs/carton



P/N	SPECIFICATION	PACKING
Ⓓ 0500-1002-02	2-pc, BLUE color. PP material	100pcs/dispenser box, 1,000pcs/carton
Ⓔ 0500-1002	2-pc, NEUTRAL color. PP material	100pcs/dispenser box, 1,000pcs/carton
Ⓕ 0500-2003	3-pc, NEUTRAL color. PP material	100pcs/dispenser box, 500pcs/carton
Ⓖ 0500-1005	5-pc, NEUTRAL color. PP material	100pcs/dispenser box, 1,000pcs/carton
Ⓗ 0500-1305	5-pc, NEUTRAL color. PP material	100pcs/dispenser box, 1,000pcs/carton



P/N	SPECIFICATION	PACKING
Ⓘ 0500-1105	5-pc, NEUTRAL color. PP material	80pcs/dispenser box, 800pcs/carton
Ⓣ 0500-1205	5-pc, NEUTRAL color. PP material	80pcs/dispenser box, 800pcs/carton




P/N	SPECIFICATION	PACKING
Ⓚ 0500-5020-16	20-pc, WHITE color. PS material	1pc/polybag, 50pcs/carton
Ⓛ 0500-5020-08	20-pc, GREEN color. PS material	1pc/polybag, 50pcs/carton
Ⓜ 0500-5020-12	20-pc, PINK color. PS material	1pc/polybag, 50pcs/carton
Ⓨ 0500-5020-02	20-pc, BLUE color. PS material	1pc/polybag, 50pcs/carton



Slide Storage Box

Designed as a convenient solution for temporary storage and transferring of microscope slides
Accommodate both the CITOTEST standard microscope slides in dimensions of 75.0mm×25.0mm×(0.8-1.2)mm and the microscope slides of 76.0mm×26.0mm×(0.8-1.2)mm in the market.

- Made of ABS material, reusable
- Snap closures/metal lock and hinged-lids provide safe operation
- A water-proof index card mounted inside the lids for simple specimen identification
- Cord sheet gaskets prevent the slides against movement shakes



P/N	COLOR	SPECIFICATION	PACKING
0500-3012-07	● Grey	12-pc,with snap closure	1pc/polybag,240pcs/carton
0500-3012-10	● Orange	12-pc,with snap closure	1pc/polybag,240pcs/carton
0500-3012-02	● Blue	12-pc,with snap closure	1pc/polybag,240pcs/carton
0500-3012-13	● Red	12-pc,with snap closure	1pc/polybag,240pcs/carton
0500-3012-16	○ White	12-pc,with snap closure	1pc/polybag,240pcs/carton
0500-3013-07	● Metallic Grey	12-pc,with snap closure	1pc/polybag,240pcs/carton
0500-3025-07	● Grey	25-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3025-10	● Orange	25-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3025-02	● Blue	25-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3025-13	● Red	25-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3025-16	○ White	25-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3026-07	● Metallic Grey	25-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3050-07	● Grey	50-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3050-10	● Orange	50-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3050-02	● Blue	50-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3050-13	● Red	50-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3050-16	○ White	50-pc,with snap closure	1pc/polybag,120pcs/carton
0500-3051-07	● Metallic Grey	50-pc,with snap closure	1pc/polybag,120pcs/carton

P/N	COLOR	SPECIFICATION	PACKING
0500-3100-07	● Grey	100-pc,with snap closure	1pc/polybag, 50pcs/carton
0500-3100-10	● Orange	100-pc,with snap closure	1pc/polybag, 50pcs/carton
0500-3100-02	● Blue	100-pc,with snap closure	1pc/polybag, 50pcs/carton
0500-3100-13	● Red	100-pc,with snap closure	1pc/polybag, 50pcs/carton
0500-3100-16	○ White	100-pc,with snap closure	1pc/polybag, 50pcs/carton
0500-3101-07	● Metallic Grey	100-pc,with snap closure	1pc/polybag, 50pcs/carton
0500-4100-07	● Grey	100-pc, w/metal lock	1pc/polybag, 50pcs/carton
0500-4100-10	● Orange	100-pc, w/metal lock	1pc/polybag, 50pcs/carton
0500-4100-02	● Blue	100-pc, w/metal lock	1pc/polybag, 50pcs/carton
0500-4100-13	● Red	100-pc, w/metal lock	1pc/polybag, 50pcs/carton
0500-4100-16	○ White	100-pc, w/metal lock	1pc/polybag, 50pcs/carton
0500-4101-07	● Metallic Grey	100-pc, w/metal lock	1pc/polybag, 50pcs/carton

Storage Cabinets

Specimen filing and storing is one of the biggest work for the pathology department. Directed by user demands and preference, CITOTEST have been manufacturing and supplying several models of storage cabinets since many years ago. The storage cabinets are constructed by Iron material or heavy-duty ABS material, in Table-Top design for less sampling and convenient transport or in the traditional standard design for mass specimen filing& storing.

BASPATH™ & BASPATH-N™ Storage Cabinets

Stackable and lockable with the series storage cabinets

Storage cabinet lid
0510-0007 /0510-0024

Storage cabinet base
0510-0006/0510-0023

Block storage cabinet
0510-0009/0510-0022

Slide storage cabinet
0510-0008/0510-0021

Each drawer accepts approx. 65 paraffin blocks

Each drawer accepts approx. 400 slides

Maglev Slide/Block Stopper
Supplied as an optional part for all the models of the CITOTEST storage cabinets, for keeping the microscope slides and paraffin blocks in upright position in the drawers.
Made of stainless steel material, 10pcs/polybag, 100pcs/box

P/N	PACKING
0510-0005	Maglev Slide/Block Stopper

LUXPATH™ Storage Cabinets

Especially designed for mass slides storage, the range up to 45, 000 slides according to each storage cabinet unit accommodates approx.4, 500 slides

Accept both the CITOTEST standard microscope slides in dimensions of 75.0mm×25.0mm×(0.8-1.2)mm and the microscope slides of 76.0mm×26.0mm×(0.8-1.2)mm in the market

- Made of heavy-duty Iron material, W/plastic powder coating, chemical resistant
- Each unit contains six interchangeable drawers, fitted with cardholder and brass knob
- Drawers slide smoothly and quietly, available with or without dividers
- Stackable and lockable storage cabinet units with the interlocking system for safe storage
- Available in white, green and tan color

The interlocking system makes it possible to stack the LuxPATH™ storage cabinets for safely storing more slides

The strong built-in separating wall between every two drawers helps realize smooth sliding.

The pre-assembled brass knob provides with convenient push-in and pull-out

Preserved grooves in front of each drawer slot accommodate long I.D labels that are supplied with the unit but not pre-installed.



COLOR	SLIDE STORAGE CABINET WITH DRAWERS DIVIDED	SLIDE STORAGE CABINET WITH DRAWERS UN-DIVIDED	STORAGE CABINET BASE WITHOUT LEG	STORAGE CABINET BASE WITH LEG
Green	0510-0037-08	0510-0031-08	0510-0038-08	0510-0032-08
White	0510-0037-16	0510-0031-16	0510-0038-16	0510-0032-16
Tan	0510-0037-15	0510-0031-15	0510-0038-15	0510-0032-15

Table-Top Storage Cabinets

Designed for conveniently storing and transferring specimens by your hand-side at your working table, instead of any other specimen holder or specimen mailer

- Accept all the regularly-sized paraffin blocks
- External Dim. L299mm×W252mm×H110mm



Mini™ Block Storage Cabinet

- Made of Iron material , with while plastic powder coating
- External Dim. L299mm×W252mm×H110mm
- Each cabinet contains 2 ABS drawers, with or without dividers
- White and orange color drawers optional

Packing: 1unit/box, 5units/carton

SPECIFICATION	BLOCK STORAGE CABINET WITH DRAWERS IN WHITE COLOR	BLOCK STORAGE CABINET WITH DRAWERS IN ORANGE COLOR
Cap. Approx.630 cassettes. W/drawers constructed w/out dividers	0510-0001-16	0510-0002-10
Cap. Approx.250 cassettes. W/drawers constructed with dividers	0510-0003-16	0510-0004-10

Slide Dispenser

Used to store the microscope slides samples

- Suitable for 25.0mm×75.0mm Microscope slides
- Constructed by white ABS material, with transparent cover in PS material



P/N	DESCRIPTION	PACKING
0500-4072	Hold up to 72 pieces of microscope slides	1pc/inner box , 20pcs/carton

MiniPLUS™ Block Storage Cabinet

- Molded from heavy-duty ABS material
- External Dim. L442mm×W250mm×H409mm
- Enclosed hand-tailored handle with polished finishing for safe as well as convenient cabinet sliding
- Each cabinet unit contains 6 drawers for up to 2,100 blocks, and marked with a plastic-covered label for identification
- All the units are stackable with built-in locking system for safely storing more blocks.



SPECIFICATION	P/N	PACKING
MiniPLUS™ block storage cabinet	0510-0050	1unit/box

Secure Archival Storage Solution NEW

Specially designed as a durable and functional product to store tissue blocks and microscope slides, dedicating to deliver the best storage solution

- Dividers included
- Unassembled, easy to set up
- Dimensions: L300mm×W295mm×H98mm
- With corrugated cardboard construction



P/N	SPECIFICATION	PACKING
0511-0003	8 drawers, holds approx 500 tissue blocks	15pcs/case
0511-0004	4 drawers, holds approx 2,000 slides	20pcs/case