



Healthcare



Healthcare Product Highlights

Innovative healthcare PPE solutions



PPE for a Sustainable Future



Made in Europe

European Sustainable & Resilient Manufacturing

Protecting our Critical Workers with Innovation & Technology



With over 25 years of heritage, Globus Group is one of the world's leading manufacturers of innovative protection solutions for healthcare and industry. Our expert team continues to work closely with healthcare providers and frontline staff. With an international footprint and manufacturing facilities based in the UK, Germany and Italy, we are providing resilient, sustainable production across Europe.



Made In Europe

Prior to the COVID-19 crisis, Globus Group had already committed to manufacturing 75% of its healthcare products in the UK and Europe. The increased demand for PPE as a result of the pandemic meant that we have been able to accelerate this strategy and significantly increase our European production, reducing our global carbon footprint. The strong and resilient networks and supply chains we have developed since launching these facilities ensure rapid delivery to our European customers.



PPE For A Sustainable Future

Our core values include a deep commitment to the long-term welfare of future generations. That means we pledge to ensure our business and operational processes are as economically, environmentally and socially sustainable as possible. Through this, we are reducing waste and energy use, protecting natural resources, supporting our communities. As well as reducing our carbon footprint, we are leading the way by investing in innovative production processes that are designed to reduce the use of single-use plastics as well as repurposing and utilising waste material.

At Globus Group, we are committed to operating as a “Greening” Company, going further to become a positive contributor to the Paris Agreement goals.

Our objectives are to:





- ✓ **Slash our carbon footprint**
- ✓ **Build a circular economy**
- ✓ **Design carbon neutral products**

European Sustainable & Resilient Manufacturing



5 European Factories Across Germany, Italy and the UK

Contents

	Respiratory Protection.....	4–12
	Hand Protection.....	13–20
	Protective Clothing.....	21–23
	Eye and Face Protection	24–38



HAIKA®

MX IIR T



Made in Europe

Tie-On Medical Face Mask Type IIR



The Haika medical face masks offer excellent protection to limit the transmission of infection.

These 3-ply masks consist of three layers of non-woven fabric including a waterproof outer layer, melt-blown middle layer and soft touch inner layer.

Key Features

- ✓ Excellent bacterial filtration efficiency 98%
- ✓ Splash-resistant
- ✓ Easy to fit and adjusted over the nose and mouth
- ✓ Non-ferrous nose wire
- ✓ Low breathing resistance and excellent comfort to wearer
- ✓ Tested in accordance of EN14683:2019



  MEDICAL DEVICE CLASS 1	Item number	HAK8000 (ties)
	Packaging	50 masks per box 40 boxes per case
	Certification	Medical Directive 93/42/ EEC and 2007/47/EC

HAIKA®

MX IIR EL



Made in Europe

Ear Loop Medical Face Mask Type IIR


The Haika medical face masks offer excellent protection to limit the transmission of infection.

These 3-ply masks consists of three layers of non-woven fabric including a waterproof outer layer, meltblown middle layer and soft touch inner layer.

Key Features

- ✓ Excellent bacterial filtration efficiency 98%
- ✓ Splash-resistant
- ✓ Easy to fit and adjusted over the nose and mouth
- ✓ Non-ferrous nose wire
- ✓ Low breathing resistance and excellent comfort to wearer
- ✓ Tested in accordance of EN14683:2019



  MEDICAL DEVICE CLASS 1	Item number	HAK8001 (earloop)
	Packaging	50 masks per box 40 boxes per case
	Certification	Medical Directive 93/42/ EEC and 2007/47/EC



2500 SERIES



4 POINT HEADSTRAP ADJUSTMENT



EXCEPTIONAL FIT



PREFORMED NOSE BRIDGE

The 2500 SERIES disposable respirators with soft inner seal offer an outstanding level of protection and fit. Suitable for protection against a broad range of particulate hazards.



2520

FFP2 Non-Valved Respirator

Key Features

- ✓ Lightweight, cup-shaped respirator
- ✓ Features a classic nose clip for a secure fit
- ✓ Soft inner nose seal
- ✓ Continuous loop head strap ensures easy-fitting

	Item number	MM25/ALP/2520
	Packaging	20 masks per box 12 boxes per case
	Certification	EN149:2001+A1:2009 FFP2 NR PPE Regulation EU 2016/425



2520V

FFP2 Valved Respirator

Key Features

- ✓ Lightweight, cup-shaped respirator
- ✓ Low-profile exhalation valve to reduce heat and moisture build up
- ✓ Features a classic nose clip for a secure fit
- ✓ Soft inner nose seal
- ✓ Continuous loop head strap ensures easy-fitting

	Item number	MM25/ALP/2520V
	Packaging	20 masks per box 12 boxes per case
	Certification	EN149:2001+A1:2009 FFP2 NR PPE Regulation EU 2016/425





ALPHA
S O L W A Y

2500 SERIES **2530V**

FFP3 Respiratory Mask With Valve

Key Features

- ✓ Lightweight, cup-shaped respirator
- ✓ Low profile exhalation valve to reduce heat and moisture buildup
- ✓ Four point head strap adjustment
- ✓ Preformed nose bridge with additional nose clip for a secure fit
- ✓ Features a full foam inner face seal for comfortable secure fit



CAT III 2797	Item number	MM25/ALP/2530V
	Packaging	10 masks per box 12 boxes per case
	Certification	EN149:2001+A1:2009 FFP2 NR PPE Regulation EU 2016/425

ALPHA
S O L W A Y

3000 SERIES **3030V**

FFP3 Disposable, Pre-formed Cup Shape Respirator

The Alpha Solway 3030V is a preformed cup shape mask offering outstanding protection coupled with exceptional levels of fit and comfort.

Key Features

- ✓ Latex, graphene and metal free
- ✓ Lightweight, cup-shaped respirator
- ✓ Four point head strap adjustment
- ✓ Electrostatic filter media
- ✓ Preformed nose bridge removes the need for regular adjustment



CAT III 2849	Item number	3030V
	Packaging	5 singles per box 20 boxes per case
	Certification	EN 149:2001+A1:2009 FFP3 NR PPE Regulation EU 2016/425



3530 SERIES

The 3500 SERIES disposable respirators offer an outstanding level of protection. Demonstrated when conducting face fit testing. Having electrostatic filter media allows low breathing resistance and a lightweight outer shell maintains excellent shape.



3530



Made in Europe

FFP3 Disposable Pre-Formed Cup Shape Respirator

Key Features

- ✓ Lightweight, cup-shaped respirator
- ✓ Four point head strap adjustment
- ✓ Preformed nose bridge removes the need for regular adjustment
- ✓ Electrostatic filter media
- ✓ Soft touch elastic
- ✓ PVC foam inner face seal



 MEDICAL DEVICE CLASS 1	Item number	3530
	Packaging	5 masks per box 20 boxes per case
	Certification	EN 149:2001+A1:2009 FFP3 NR EN 14683:2019+AC:2019 FFP3 NR + PPE Reg MDD (UK) 93/42/EEC & 2007/47/EC MDR (EU) 2017/745



3530V



Made in Europe

FFP3 Disposable Pre-Formed Cup Shape Respirator

Key Features

- ✓ Lightweight, cup-shaped respirator
- ✓ Four point head strap adjustment
- ✓ Preformed nose bridge removes the need for regular adjustment
- ✓ Electrostatic filter media
- ✓ Soft touch elastic
- ✓ PVC Foam Inner Face Seal
- ✓ Low-profile exhalation valve helps reduce heat and moisture build up



	Item number	3530V
	Packaging	5 masks per box 20 boxes per case
	Certification	EN149:2001+A1:2009 FFP3 NR PPE Regulation EU 2016/425



ALPHA
S O L W A Y

ALPHA S

FFP3 Pre-Formed Cup Shape Respirator

The Alpha Solway S is a preformed cup shape mask specifically designed to ensure an excellent fit for the smaller face.

Key Features

- ✓ Latex, graphene and metal-free
- ✓ Preformed nose bridge removes the need for regular adjustment
- ✓ Unique inner face seal helps maintain a close and snug fit around the chin area
- ✓ Four point adjustable head strap for a secure fit
- ✓ Easy to fit for first time users



	Item number	BWG069
	Packaging	5 masks per box 20 boxes per case
	Certification	EN149:2001+A1:2009 FFP3 NR (EU) 2016/425

ALPHA
S O L W A Y

ALPHASOL+ Model: AP-3




Made in Europe

FFP3 Vertical Fold Flat Respirator

The Alphasol+ is a vertical fold flat disposable respirator designed with an inner face seal to offer an excellent fit for all face shapes.

Key Features

- ✓ Lightweight, vertical, fold flat respirator
- ✓ Flexible nose bridge
- ✓ Inner seal around the nose area for an effective fit
- ✓ Electrostatic filter media
- ✓ Valve-free
- ✓ Aluminium nose clip
- ✓ Latex-free
- ✓ Premium quality straps



	Item number	MFAS/ALP/AP3
	Packaging	20 singles per box 12 boxes per case
	Certification	EN 149:2001+A1:2009 FFP3 NR EN 14683:2019+AC:2019 PPE Regulation EU 2016/425



Disposable Vertical Fold Flat Respirator

The AlphaFold+ is a vertical fold flat disposable respirator designed with an inner face seal, adjustable straps and two size options to offer an excellent fit for all face shapes.

Key Features

- ✓ Lightweight, fold flat respirator
- ✓ Available in two size options
- ✓ Inner seal around the nose area and fully adjustable straps
- ✓ Protection options: FFP2 - FFP3
- ✓ Available in smaller size (ASFS-2, ASFS-2V, ASFS-3, ASFS-3V)
- ✓ Valved and unvalved

ASF-2



ASF-3



	Item number	ASF-2/ASFS-2	ASF-2V/ASFS-2V
	Packaging	20 masks per box 12 boxes per case	15 masks per box 12 boxes per case
	Certification	EN 149:2001+A1:2009 • FFP2 NR PPE Regulation EU 2016/425	EN 149:2001+A1:2009 • FFP2V NR PPE Regulation EU 2016/425

	Item number	ASF-3/ASFS-3	ASF-3V/ASFS-3V
	Packaging	20 masks per box 12 boxes per case	15 masks per box 12 boxes per case
	Certification	EN 149:2001+A1:2009 • FFP3 NR PPE Regulation EU 2016/425	EN 149:2001+A1:2009 • FFP3V NR PPE Regulation EU 2016/425



The HX-SERIES is a horizontal fold flat disposable respirator designed to offer an excellent fit.

Key Features

- ✓ Lightweight, fold flat respirator
- ✓ Flexible nose bridge
- ✓ Electrostatic filter media
- ✓ Tapered shape for enhanced fit
- ✓ Protection options: FFP2 - FFP3
- ✓ Available in larger size (HXL-2, HXL-3)
- ✓ Latex and silicone-free



	Item number	HX-2/HXL-2	HX3/HXL3
	Packaging	20 masks per box 12 boxes per case	20 masks per box 12 boxes per case
	Certification	EN 149:2001+A1:2009 • FFP2 NR PPE Regulation EU 2016/425 Medical Directive 93/42/ EEC and 2007/47/EC EN14683:2019+AC:2019	EN 149:2001+A1:2009 • FFP3 NR PPE Regulation EU 2016/425 Medical Directive 93/42/ EEC and 2007/47/EC EN14683:2019+AC:2019



HAIKA | LX Series



Made in Europe

Disposable Fold Flat Respirator

The LX-SERIES is a horizontal fold flat disposable respirator designed to offer an excellent fit.

Key Features

- ✓ Lightweight, fold flat respirator
- ✓ Flexible nose bridge
- ✓ Electrostatic filter media
- ✓ Protection options: FFP2 - FFP3
- ✓ Premium quality straps ensure stability on the face
- ✓ Non-ferrous nose bridge



	Item number	HAK4000	HAK4002
	Packaging	50 masks per box 12 boxes per case	50 masks per box 12 boxes per case
	Certification	Medical Directive 93/42/EEC and 2007/47/EC. EN 149:2001+A1:2009 FFP2 NR	Medical Directive 93/42/EEC and 2007/47/EC. EN 149:2001+A1:2009 FFP3 NR

ALPHA SOLWAY | HXC-SERIES



Made in Europe

Disposable Fold Flat Respirator with a Clear Panel

The HXC-Series are horizontal fold flat disposable respirators offering excellent protection to limit the transmission of infection with the addition of a clear anti-fog front panel, positioned to prevent reflection and make lip reading easier.

Key Features

- ✓ Lightweight, fold flat respirator
- ✓ Flexible nose bridge
- ✓ Valve-free
- ✓ Protection options: FFP2 - FFP3
- ✓ Latex and silicone-free
- ✓ Soft touch elastic
- ✓ Non-sterile



	Item number	HXC-2	HXC-3
	Packaging	20 masks per box 12 boxes per case	20 masks per box 12 boxes per case
	Certification	Conforms to EN 149:2001+A1:2009 FFP2 NR & tested in accordance of transparent face mask technical specification: 1st April 2021 PPE Regulation EU 2016/425 Medical Directive 93/42/EEC and 2007/47/EC	Conforms to EN 149:2001+A1:2009 FFP2 NR & tested in accordance of transparent face mask technical specification: 1st April 2021 PPE Regulation EU 2016/425 Medical Directive 93/42/EEC and 2007/47/EC



HAIKA | MX C



Made in Europe


Face Mask with Clear Panel

The Haika MX C Transparent Face Mask offers excellent protection to limit the transmission of infection with the addition of a clear anti-fog front panel, positioned to prevent reflection and make lip reading easier.

Key Features

- ✓ Excellent bacterial filtration efficiency 98%
- ✓ Splash-resistant
- ✓ Valve-free
- ✓ Latex and silicone-free
- ✓ Hypoallergenic



 MEDICAL DEVICE CLASS 1	Item number	HAK6000
	Packaging	50 masks per box 12 boxes per case
	Certification	Medical Directive 93/42/EEC and 2007/47/EC. Tested in accordance of Transparent Face Mask Technical specification: 1st April 2021

Benchmark | BM-Series

The Benchmark BM-Series cup-shaped respirators offer outstanding protection coupled with exceptional levels of fit and comfort. The elasticated straps are colour-coded for ease of identification.


BM-2

FFP2 Respiratory Mask Non-Valved

The BM-2 non-reusable respirator features a soft foam nose seal, blue twin elastic strap that are welded and not stapled for strength. Colour-coded for ease of identification.

Key Features

- ✓ Soft inner nose seal
- ✓ Blue colour-coded noseband and strap
- ✓ Low breathing resistance
- ✓ Staple-free
- ✓ Soft inner lining

 CAT III CE 0086	Item number	BMR00020DZ
	Packaging	20 masks per box 400 per case
	Certification	EN149: 2001 + A1: 2009 FFP2 NR PPE Regulation EU 2016/425





Face Fit Testing

COMPULSORY TRAINING

Training in the proper use and care of respirators has been mandatory in the COSHH (Control of Substances Hazardous to Health) regulations since 1989. The 2002 revision of the COSHH regulations and the Approved Code of Practice (ACoP) now require that fit testing of Respiratory Protective Equipment (RPE) be carried out. Only respirators that are loose fitting and rely on airflow through the device are exempt e.g. powered helmets.

LEADING PROVIDER

Alpha Solway has staff that were the first to obtain the 'Fit2Fit' accreditation. 'Fit2Fit' is recognised within INDG479 Fit Testing Guidance Notes.

The scheme confirms the competency of a particular individual to perform face piece fit testing, by ensuring they have passed an industry recognised exam and demonstrated a thorough knowledge of how to fit test in one or more of the test methods.

Face Fit Testing Kit

The Alpha Solway Qualitative Face Fit Testing Kits enable wearers of negative pressure half masks and disposable respirators to be Fit Tested correctly and economically. Four variations of the Kit are available depending on whether bitter or sweet sensitivity solutions are required and with a choice of cases.

All are easy to use and supplied with full instructions. Spare parts are available to order upon request.



QLTK1	QLTK1P	QLTK1S	QLTK1SP
1 x Hood	1 x Hood	1 x Hood	1 x Hood
12 x 2.5ml Bitter Sensitivity Solutions	12 x 2.5ml Bitter Sensitivity Solutions	12 x 2.5ml Sweet Sensitivity Solutions	12 x 2.5ml Sweet Sensitivity Solutions
12 x 2.5ml Bitter Test Solutions	12 x 2.5ml Bitter Test Solutions	12 x 2.5ml Sweet Test Solutions	12 x 2.5ml Sweet Test Solutions
2 x Nebulisers	2 x Nebulisers	2 x Nebulisers	2 x Nebulisers
Instructions	Instructions	Instructions	Instructions
Supplied in a Cardboard Case	Supplied in a Durable Plastic Case	Supplied in a Cardboard Case	Supplied in a Durable Plastic Case



Train The Trainer

Alpha Solway Fit2Fit Accredited experts can also conduct 'Train the Trainer' courses. These are designed to enable companies to carry out the Fit Testing of wearers with complete confidence. Contact us for more details.



HAIKA | NX510

Single-Use, Powder-Free, Nitrile Examination Glove

Lightweight, ambidextrous, single-use medical nitrile examination gloves coloured light blue.

Key Features

- ✓ A textured surface allows for a secure grip
- ✓ Supreme comfort and sensitivity
- ✓ Effective barrier protection
- ✓ Beaded cuff
- ✓ Powder-free and latex-free

						Packaging: 100 gloves per dispenser 10 dispensers per case
--	--	--	--	--	--	---

Medical Directive 93/42/EEC and 2007/47/EC PPE Regulation EU 2016/425, MDR 2017/745 Class I
EN ISO 374-1: 2016 Type B (KPT) /EN ISO 374-5: 2016 (Virus)/ EN455 (Parts 1 to 4), EN421:2010
(Excluding clause 4.3)

NX510

SIZE	CODE
XS	HAK5100
S	HAK5101
M	HAK5102
L	HAK5103
XL	HAK5104
XXL	HAK5105

NX520

SIZE	CODE
XS	HAK5200
S	HAK5201
M	HAK5202
L	HAK5203
XL	HAK5204
XXL	HAK5205

Also available
NX520
200 gloves per dispenser
XL and XXL 180 single gloves
10 boxes per case



Benchmark | BMG462

Single-Use, Powder-Free, Nitrile Glove

Lightweight, powder-free, single-use gloves with textured fingertips. Provides a very good alternative to vinyl gloves.

Key Features

- ✓ Textured finish on the fingertips to enhance grip
- ✓ Comfortable with a 'second skin' feel
- ✓ Powder-free and latex-free
- ✓ Approved for food handling
- ✓ Beaded cuff

					Packaging: 200 gloves per dispenser 10 dispensers per case
--	--	--	--	--	---

PPE Regulation EU 2016/425 EN ISO 374-1: 2016 Type B (KTP)/EN ISO 374-5: 2016
(Virus)/EN455 (Parts 1 to 4) Medical Directive 93/42/EEC MDR 2017/745 Class I

BMG462

SIZE	CODE
6½-7	BMG00462DD
7½-8	BMG00462DF
8½-9	BMG00462DH
9½-10	BMG00462DJ
10½-11	BMG00462DL

BMG460

SIZE	CODE
6½-7	BMG00460DD
7½-8	BMG00460DF
8½-9	BMG00460DH
9½-10	BMG00460DJ
10½-11	BMG00460DL

Also available
BMG460
100 gloves per dispenser
10 boxes per case





Benchmark TM

BMG450

Nitrile, Powder-Free, Single-Use Protection

High performance, powder-free, nitrile, single-use disposable glove. Significantly stronger alternative to latex.

Key Features

- ✓ Significantly stronger alternative to latex
- ✓ Comfortable with a 'second skin' feel
- ✓ Resistant to oils, acids and many chemicals
- ✓ Textured finish on palm for improved grip and touch
- ✓ Powder-free and latex-free

CAT III CE 0120	EN ISO 374-1: 2016 TYPE B JKT	EN ISO 374-5: 2016 VIRUS	Packaging: 100 gloves per box 10 dispensers per box
------------------------------	----------------------------------	-----------------------------	--

SIZE	CODE
6½-7	BMG00450DD
7½-8	BMG00450DF
8½-9	BMG00450DH
9½-10	BMG00450DJ

PPE Regulation EU 2016/425 EN ISO 374-1: 2016 Type B (JKT)/ EN ISO 374-5: 2016 / EN455 (Parts 1 to 4)



SKYTEC TM

Teal TM

Nitrile, Powder-Free, Single-Use Protection

Premium performance, powder-free, nitrile, single-use disposable glove. Extremely strong and offers greater chemical and puncture resistance than regular nitrile disposables.

Key Features

- ✓ 100% premium quality nitrile, thickness 0.12mm
- ✓ Resistant to oils, acids and many chemicals
- ✓ Up to four times the puncture resistance of latex
- ✓ Powder-free, non-sterile
- ✓ Lightweight and comfortable with a 'second skin' feel

CAT III CE 0598	EN455 MD MEDICAL DEVICE	EN ISO 374-1: 2016 TYPE B JKT	EN ISO 374-5: 2016 VIRUS	Packaging: 100 singles per box 10 boxes per case
------------------------------	--------------------------------------	----------------------------------	-----------------------------	---

SIZE	CODE
6/XS	SKYTEAL0
7/S	SKYTEAL1
8/M	SKYTEAL2
9/L	SKYTEAL3
10/XL	SKYTEAL4
11/XXL	SKYTEAL5

Medical Directive 93/42/EEC & 2007/47/EC MDR 2017/745 Class I PPE Regulation EU 2016/425 EN ISO 374-1: 2016 Type B (JKPT)/EN ISO 374-5: 2016 (Virus)/ EN455 (Parts 1 to 4)



FOOD SAFE





Iris™



Nitrile, Powder-Free, Single-Use Examination Glove

High performance, powder-free, nitrile, single-use examination glove. Significantly stronger alternative to latex and food safe.

Key Features

- ✓ 100% premium quality nitrile, thickness 0.08mm
- ✓ Powder-free
- ✓ Comfortable with a 'second skin' feel
- ✓ Resistant to oils, acids and many chemicals
- ✓ Rolled cuff for secure fit
- ✓ Non-sterile
- ✓ Chemotherapy approved
- ✓ Conforms to EN455 (all parts), ASTM D6319-10 and D6124-06

CAT III 2777	EN ISO 374-1: 2016 TYPE B 	EN ISO 374-5: 2016 		Packaging: 100 singles per box 10 boxes per case
---------------------	-------------------------------	------------------------	--	---

SIZE	CODE
7/S	SKYIRIS1
8/M	SKYIRIS2
9/L	SKYIRIS3
10/XL	SKYIRIS4
11/XXL	SKYIRIS5

Medical Directive 93/42/ EEC and 2007/47/EC PPE MDR 2017/745 Class I Regulation EU 2016/425 EN ISO 374-1: 2016 Type B (KTP) /EN ISO 374-5: 2016 (Virus) / EN455 (Parts 1 to 4) Conforms to EN455 (all parts), ASTM D6319-10 and D6124-06



Utah™



Nitrile, Powder-Free, Single-Use Examination Glove

High performance, powder-free, nitrile, single-use examination glove. Significantly stronger alternative to latex and food safe.

Key Features

- ✓ Significantly stronger alternative to latex
- ✓ 100% premium quality nitrile, thickness 0.10mm
- ✓ Comfortable with a 'second skin' feel
- ✓ Resistant to oils, acids and many chemicals
- ✓ Textured finish on fingertips for improved grip and touch
- ✓ Powder-free, non-sterile
- ✓ Conforms to EN455 (all parts), ASTM D6319-10 and D6124-06

CAT III 2777	 MEDICAL DEVICE CLASS 1	EN ISO 374-1: 2016 TYPE B 	EN ISO 374-5: 2016 		Packaging: 100 singles per box 10 boxes per case
---------------------	-------------------------------	-------------------------------	------------------------	--	---

SIZE	CODE
7/S	SKYUTAH1
8/M	SKYUTAH2
9/L	SKYUTAH3
10/XL	SKYUTAH4
11/XXL	SKYUTAH5

Medical Directive 93/42/EEC and 2007/47/EC MDR 2017/745 Class I PPE Regulation EU 2016/425 EN ISO 374-1: 2016 Type B (KTP) /EN ISO 374-5: 2016 (Virus) / EN455 (Parts 1 to 4) ASTM D6319-10 and D6124-06

i-con™
Glove System
Part of the
ICON range





I-CON™ is a colour-coded glove system that provides enhanced chemical protection for workers in the cleaning and hygiene sector.

As part of the range, it includes four nitrile single-use gloves. The different colour of the gloves represent the four designated aspects of cleaning and hygiene as per the British Institute of Cleaning Science (BICSc) guidelines for helping to eliminate cross contamination within end-user sites. Whether used for a cleaning application as part of the Globus I-CON system or for light chemical protection application within other industries, these four high performance Skytec gloves are both durable and comfortable.




		SKYTEC™ RHODE	SKYTEC™ TEXAS	SKYTEC™ ALASKA	SKYTEC™ UTAH
ICON™ USE		Washroom Areas	Clinical, Treatment & Isolation Areas	Kitchen Areas	General Areas (Excluding washrooms, clinical areas and kitchens)
XS	SKY980	SKY1000	SKY990	-	
S	SKY981	SKY1001	SKY991	SKYUTAH1	
M	SKY982	SKY1002	SKY992	SKYUTAH2	
L	SKY983	SKY1003	SKY993	SKYUTAH3	
XL	SKY984	SKY1004	SKY994	SKYUTAH4	
XXL	-	-	-	SKYUTAH5	
10 BOXES OF 100					



Key Features

- ✓ Excellent fingertip sensitivity
- ✓ Comfortable with a 'second skin' feel
- ✓ Textured finish on fingertips for improved grip
- ✓ Powder-free, non-sterile
- ✓ AQL 1.5
- ✓ CE category III approved
- ✓ Approved for food handling
- ✓ Latex free - ideal for people with Type I latex sensitivity
- ✓ Conforms to ASTM D6319-10 and D6124-06
- ✓ Length 240mm



SUITABLE APPLICATIONS

The Globus I-CON™ system is suitable for using as part of a managed colour-coded cleaning programme within a wide number of facilities including:



Manufactured from nitrile, which offers additional strength over latex or vinyl, these chemical gauntlets are longer and offer more resistance to chemicals than shorter, disposable gloves.




Key Features

- ✓ Flock-lined for extra comfort and to ease donning
- ✓ Textured palm and fingertips to enhance grip in wet conditions
- ✓ CE category III approved
- ✓ Chlorinated and sanitised to reduce odours
- ✓ Straight cuff
- ✓ Latex-free - ideal for people with Type I latex sensitivity
- ✓ Thickness: 0.38mm
- ✓ Length: 330mm

SUITABLE APPLICATIONS



ICON™ USE	SKYTEC™ DAKOTA	SKYTEC™ FRISCO
	Kitchen Areas	General Areas (Excluding washrooms, clinical areas and kitchens)
XS	SKY960	SKY950
S	SKY961	SKY951
M	SKY962	SKY952
L	SKY963	SKY953
XL	SKY964	SKY954
CASE QTY – 120		



For further information please visit www.globus.co.uk/icon





“ HAIKA provides a quality range of gloves. In particular, it has produced a synthetic glove with the handling characteristics of the best latex gloves. It has become my glove of choice in clinical practice. ”

MEMBER AND FORMER PRESIDENT OF THE COUNCIL OF THE BRITISH ASSOCIATION OF AESTHETIC PLASTIC SURGEONS

HAIKA | Intra


Sterile

Sterile, Latex, Powder-Free Surgical Undergloves

Designed to reduce risks when performing orthopaedic surgery, HAIKA™ INTRA is a soft and elastic dark green coloured latex underglove. INTRA can be paired with any surgical overglove for enhanced protection during intensive surgical procedures.

Key Features

- ✓ Enhanced safety
- ✓ Reduced hand fatigue
- ✓ Excellent double donning
- ✓ Powder-free



EN455 1,2 & 3
Medical Directive 93/42/EEC and 2007/47/EC

SIZE	CODE
5½	HAK6055
6	HAK6060
6½	HAK6065
7	HAK6070
7½	HAK6075
8	HAK6080
8½	HAK6085
9	HAK6090





Callisto

Sterile

Sterile, Latex, Powder-Free Surgical Gloves

HAIKA™ CALLISTO is a natural coloured, powder-free surgical glove with enhanced puncture resistant properties. Expertly manufactured from natural rubber latex using an advanced deproteinised process, CALLISTO features unequivocal elasticity and contains significantly reduced levels of proteins.

Key Features

- ✓ Secure control
- ✓ Supreme donning and comfort
- ✓ Reduced hand fatigue
- ✓ Powder-free
- ✓ Chemical resistant



CAT III

 2797

MEDICAL DEVICE
 CLASS 2A

Packaging:
 50 pairs per box
 4 boxes per case

EN455 1,2 & 3
Medical Directive 93/42/EEC and 2007/47/EC

SIZE	CODE
5½	HAK2055
6	HAK2060
6½	HAK2065
7	HAK2070
7½	HAK2075
8	HAK2080
8½	HAK2085
9	HAK2090



Syntis

Sterile

Sterile, Synthetic, Non-Latex, Powder-Free Surgical Gloves

HAIKA™ SYNTIS is a natural coloured, durable powder-free surgical glove delivering an incredible feel that far exceeds the comfort experience of a regular synthetic surgical glove. The special soft synthetic formulation contributes to an incredibly efficient donning experience for the user, even when double donning and damp hand donning.

Key Features

- ✓ Excellent fingertip sensitivity with consistent grip
- ✓ Syntis features a multi-layered polymer coating
- ✓ Expertly manufactured and anatomically designed
- ✓ Chemotherapy approved
- ✓ Powder-free



CAT III

 2797

MEDICAL DEVICE
 CLASS 2A

Packaging:
 50 pairs per box
 4 boxes per case

EN455 1,2 & 3
Medical Directive 93/42/EEC and 2007/47/EC

SIZE	CODE
5½	HAK4055
6	HAK4060
6½	HAK4065
7	HAK4070
7½	HAK4075
8	HAK4080
8½	HAK4085
9	HAK4090



HAIKA | Syntis PX

Sterile

Sterile, Polychloroprene Surgical Gloves

HAIKA™ SYNTIS PX is a latex-free microtextured surgical glove. Extremely fine, yet strong, it's also an exceptionally tactile glove that provides a remarkable fit and feel. The textured surface enhances slip resistance and the synthetic formulation allows for easy donning and doffing.

Key Features

- ✓ Multi-layered polymer coating which allows donning when double donning
- ✓ Excellent fingertip sensitivity with consistent grip
- ✓ Expertly manufactured and anatomically designed
- ✓ Chemotherapy approved
- ✓ Powder-free



CAT III
 2797
 MEDICAL DEVICE CLASS 2A
Packaging:
 50 pairs per box
 4 boxes per case

SIZE	CODE
5½	HAK7055
6	HAK7060
6½	HAK7065
7	HAK7070
7½	HAK7075
8	HAK7080
8½	HAK7085
9	HAK7090

EN455 1,2 & 3
 Medical Directive 93/42/EEC and 2007/47/EC

HAIKA | PURO | CX335™

Cleanroom

Nitrile, Powder-Free, Non-Sterile, Ambidextrous Gloves

ISO Class 5 cleanroom gloves, 100% nitrile and powder-free.

Key Features

- ✓ Beaded cuff
- ✓ High puncture resistance
- ✓ Non-sterile
- ✓ Ambidextrous
- ✓ CX334: ISO CLASS 5 compatible (blue colour)
- ✓ CX335: ISO CLASS 5 compatible (white colour)
- ✓ 300mm long by 0.13mm thick



CAT III
 0120
 MEDICAL DEVICE
 EN ISO 374-1: 2016 TYPE B
 EN ISO 374-5: 2016
 JKT
Packaging:
 100 pairs per box
 10 boxes per case

SIZE	CODE (CX334)	CODE (CX335)
6/XS	-	PUR3350
7/S	PUR3341	PUR3351
8/M	PUR3342	PUR3352
9/L	PUR3343	PUR3353
10/XL	PUR3344	PUR3354

Medical Directive 93/42/EEC & 2007/47/EC EN ISO 374-1: 2016
 Type B (JKT) / EN ISO 374-5: 2016 / EN455 (Parts 1 to 4)



HAIKA | HKC01



Made in Europe

Reusable Surgical Gown

Key Features

- ✓ Ergonomic fit for freedom of movement
- ✓ Secure back neck studded adjustment
- ✓ Elasticated soft touch ribbed cuffs
- ✓ Thumb loops for added security
- ✓ One size fits all wraparound design for enhanced coverage
- ✓ Machine washable up to 75°C for 100 washes 'under testing'



Rear view



		Packaging: 1 Single per bag
MEDICAL DEVICE CLASS 1		

ITEM NUMBER
HKC0001DZ

EN13795-1
Medical Directive 93/42/EEC and 2007/47/EC

HAIKA | HKC02



Made in Europe

Reusable Surgical Gown

Key Features

- ✓ Ergonomic fit for freedom of movement
- ✓ Secure back neck studded adjustment
- ✓ Elasticated soft touch ribbed cuffs
- ✓ Thumb loops for added security
- ✓ One size fits all wraparound design for enhanced coverage
- ✓ Machine washable up to 73°C for 100 washes
- ✓ Fluid repellent PVC coating with soft comfortable feel



Rear view



		Packaging: 1 Single per bag
MEDICAL DEVICE CLASS 1		

ITEM NUMBER
HKC0002EZ

EN13795-1
Medical Directive 93/42/EEC and 2007/47/EC



Limited Life Coverall Type 3/4

Alphachem X150B is a light, durable and cost effective disposable chemical protective garment manufactured using lightweight HDPE fabric. It offers chemical resistance against splashes and sprays of hazardous chemicals in Type 3 and 4 applications.

Key Features

- ✓ Designed to provide increased wearer comfort, exceptional levels of safety and protection, with superior durability
- ✓ Elasticated wrists, ankles and waist 2 way front zipper with resealable extended storm flap
- ✓ Integrated thumb loops
- ✓ Taped internal overlapped seams provide a strong and effective chemical barrier



CAT III CE 0493	EN14605: 2005 +A1:2009 00000X	EN14605 +A1:2009 Type 3	EN14605 +A1:2009 Type 4	ENISO 13982-1:2004 +A1:2009 Type 5	EN ISO 374-5:2016 Type 3: 4-8; 5-8; 6-8	EN511: 2016 TIL Class 1	EN1149-1	Packaging: 40 per case
-------------------------------------	----------------------------------	----------------------------	----------------------------	---------------------------------------	--	----------------------------	----------	----------------------------------

Conform type 3/4 PPE Regulation EU 2016/425

SIZE	CODE
Yellow Small	CDX1B/HX150B/YE1 X150-B-005
Yellow Medium	CDX1B/HX150B/YE2 X150-B-005
Yellow Large	CDX1B/HX150B/YE3 X150-B-005
Yellow X/Large	CDX1B/HX150B/YE4 X150-B-005
Yellow XX/Large	CDX1B/HX150B/YE5 X150-B-005
Yellow XXX/Large	CDX1B/HX150B/YE6 X150-B-005

Limited Life Chemical Coverall

The Alphachem X350 is constructed from a unique non-woven composite, engineered to offer a high level of barrier to the permeation of organic, inorganic chemicals and infective agents.

Key Features

- ✓ Manufactured using polyethylene and EVOH barrier composite, bonded structural spunbonded and polypropylene base layer suitable for higher risk applications
- ✓ Double zip system and grab tag to support with an emergency. Easy to hold when wearing gloves
- ✓ Double sleeve system with inner thumb loop ensuring arm protection to support when arm is raised



CAT III CE 0493	EN13034: 2005 +A1:2009 00000X	EN14605 +A1:2009 Type 3	EN14605 +A1:2009 Type 4	ENISO 13982-1:2004 +A1:2010 Type 5	EN ISO 374-5:2016 Type 3: 4-8; 5-8; 6-8	EN1073-2: 2003 TIL Class 1	EN1149-5: 2009	Packaging: 40 per case
-------------------------------------	----------------------------------	----------------------------	----------------------------	---------------------------------------	--	-------------------------------	----------------	----------------------------------

Conform type 3/4 PPE Regulation EU 2016/425

SIZE	CODE
X350 Orange Medium	ASC00002GF
X350 Orange Large	ASC00002GH
X350 Orange X-Large	ASC00002GJ
X350 Orange XX-Large	ASC00002GL
X350 Orange XXX-Large	ASC00002GN



Limited Life Coverall Type 5/6

Alphashield S2200 is a limited life coverall with bound seam made from polypropylene and polyethylene. It is designed to protect workers from hazardous substances or sensitive products and processes from contamination.

Key Features

- ✓ Lightweight microporous polyethylene (PE) laminate
- ✓ Bound seam for additional protection and strength
- ✓ 2 way front zipper with resealable storm flap
- ✓ Thumb loop ensures arm protection remains in place when arm is raised
- ✓ Elasticated ankles, wrists and waist to enhance fit
- ✓ Innovative, three piece elasticated hood construction provides added comfort, durability and fit



CAT III
CE
0493

EN13034: 2005 +A1:2009
Type 6

ENISO 13982-1:2004 +A1:2010
Type 5

EN14126: 2003
Type 5-B: G-B

EN1073-2: 2005
TIL Class 1

EN1149-1

Packaging:
50 per case

SIZE	CODE
White Small	PDS2B/HS2200/WH1 S2200
White Medium	PDS2B/HS2200/WH2 S2200
White Large	PDS2B/HS2200/WH3 S2200
White X/Large	PDS2B/HS2200/WH4 S2200
White XX/Large	PDS2B/HS2200/WH5 S2200
White XXX/Large	PDS2B/HS2200/WH6 S2200

Conform type 5/6 PPE Regulation EU 2016/425

Benchmark BMC-02

Limited Life Coverall SMMS Type 5/6

Benchmark BMC-02 is a limited life coverall made from lightweight spun-bonded polypropylene with an outer breathable film. It is designed to protect workers from hazardous substances or sensitive products and processes from contamination.

Key Features

- ✓ Microporous laminated material
- ✓ Elasticated waist
- ✓ Overlocked stitched seams
- ✓ Inset sleeve for ease of movement
- ✓ Available in white

CAT III
CE
0624

EN13034: 2005 +A1:2009
Type 6

ENISO 13982-1:2004 +A1:2010
Type 5

EN14126: 2003
Type 5-B: G-B

EN11: 2016
TIL Class 2

EN1149-1: 2008

Packaging:
40 per case

SIZE	CODE
M	BMC00002AF
L	BMC00002AH
XL	BMC00002AJ
XXL	BMC00002AL

Conform type 5/6 PPE Regulation EU 2016/425





Riley®

Kosma™



DIELECTRIC



TEMPLEX



SOFT TPE

Safety Glasses

A frameless design provides a sports-style look with enhanced peripheral vision. The soft TPE nose bridge and temple tips ensure enhanced comfort and a secure fit for a wide range of face sizes.



Clear



Grey

Key Features

- ✓ Anti-scratch and anti-fog dual coated lenses
- ✓ Slim wraparound style
- ✓ Soft TPE nose piece for enhanced comfort and fit
- ✓ Flexible temple tips
- ✓ 100% metal-free
- ✓ Supplied with neck cord

CAT II CE	Packaging	10 per box / 200 per case
	Certification	EN 166:2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR ○	RLY00321	2C-1.2 e 1FT CE	e EN 166 FT CE /Z87+	Impacts, UV radiations	Riley Anti-scratch and Anti-fog
GREY ●	RLY00322	5-2.5 e 1FT CE	e EN 166 FT CE /Z87+	Impacts, UV and Solar radiations	

• Product weight 26g | • Base lens 9mm



Riley® | Ligerá™



SPORTS
STYLE



ULTRA
LIGHTWEIGHT



AERO TEMPLE

Ultra-Lightweight Safety Glasses

Super lightweight Ligerá™ combines maximum wearer fit and comfort with a sleek sports design. With its high impact resistance, anti-scratch and anti-fog coatings inside and out, Ligerá™ withstands exposure to a host of typical worker environment conditions.



Clear



Blue



Grey

Key Features

- ✓ Anti-scratch and anti-fog dual coated lenses
- ✓ Highly curved lenses offer full wraparound protection and optimum vision
- ✓ High impact resistance under extremes of temperature (-5~55°C)
- ✓ Ultra-lightweight 16.5g

CAT II CE	Packaging	10 per box / 200 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425


VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR ○	RLY00490	2C-1.2 e 1FT CE Z87+U6 S/IO 0	e EN 166 FT CE / Z87+ / AS/NZS	Impacts, UV radiations	Riley Anti-Scratch / Anti-fog
GREY ●	RLY00491	5-2.5 e 1FT CE Z87+U6 S L2.5/I2	e EN 166 FT CE / Z87+ / AS/NZS	Impacts, UV and Solar radiations	
DARK GREY ●	RLY00492	5-3.1 e 1FT CE Z87+U6 S L3/I3	e EN 166 FT CE / Z87+ / AS/NZS	Impacts, UV and Solar radiations	
BLUE ●	RLY00493	2C-1.4 e 1FT CE Z87+U6 S / I1	e EN 166 FT CE / Z87+ / AS/NZS	Impacts, UV radiations	

• Product weight 16.5g | • Base lens 10mm



Safety Glasses

Stream Evo™ safety glasses provide extreme comfort and 'wearability' with a range of lens options to suit a variety of conditions.



SPORTS STYLE



DIELECTRIC



SOFT TPE



ADJUSTABLE NOSEBRIDGE



Clear

Key Features

- ✓ Single wraparound polycarbonate lens
- ✓ Soft, ventilated TempFLEX™ temple tips provide a comfortable, secure fit
- ✓ Interchangeable and removable adjustable headband
- ✓ High performance TECTON 400 anti-scratch and anti-fog coatings applied to BOTH sides of the lens
- ✓ Soft, adjustable, anti-slip nose bridge
- ✓ Slim temples enhance comfort when worn with ear defenders
- ✓ Anti-reflective and high contrast lens options
- ✓ 100% metal-free
- ✓ Supplied with a microfibre storage pouch and premium cord
- ✓ Also available in a smaller size - Stream™ Evo Small
- ✓ Also available without plastic packaging - see Stream™ (Evo) Eco

	Packaging	10 per box 200 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425



Blue

TECTON 400



Grey

TECTON 400



Amber

TECTON 400



Anti-reflective

ARIDA



Twilight

ARIDA



LED

TECTON 400

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR	○	RLY00361	2C-1.2 e 1FTKN CE / Z87+ U6 IO 0	e EN 166 FT CE / Z87+ AS / NZS 1337.1 Lic, No. 1502	Impacts, UV radiations
GREY	●	RLY00362	5-2.5 e 1FTKN CE / Z87+ U6 L2.5 I 2	e EN 166 FT CE / Z87+ AS / NZS 1337.1 Lic, No. 1502	Impacts, UV and Solar radiations
BLUE	●	RLY00365	2-1.4 e 1FTKN CE / Z87+ U6	e EN 166 FT CE / Z87+	Impacts, UV radiations
AMBER	●	RLY00371	2C-1.2 e 1FTKN CE / Z87+ U6	e EN 166 FT CE / Z87+	Impacts, UV radiations
LED	●	RLY00369	2C-1.4 e 1FTKN CE / Z87+ U6 L1.3	e EN 166 FT CE / Z87+ AS / NZS 1337.1 Lic, No. 1502	Impacts, UV radiations
ANTI-REFLECTIVE	●	RLY00363	2C-1.2 e 1FT CE / Z87+ U6 L1.7 I 1	e EN 166 FT CE / Z87+	Impacts, UV and Solar radiations
TWILIGHT	●	RLY00364	5-1.7 e 1FT CE / Z87+ U6 L1.7 I 1	e EN 166 FT CE / Z87+ AS / NZS 1337.1 Lic, No. 1502	Impacts, UV and Solar radiations
EVO ECO		RLY00370	2C-1.2 e 1FTKN CE / Z87+ U6 IO 0	e EN 166 FT CE / Z87+ AS / NZS 1337.1 Lic, No. 1502	Impacts, UV radiations
ACCESSORY PACK		RLYA0017	N/A	N/A	N/A

• Product weight 23g | • Base lens 9.5mm



Safety Glasses

Stream Evo™ Small safety glasses have all the great features of Stream Evo but offer a wide field of vision with a narrower fit for a smaller head size/shapes.



SPORTS STYLE



DIELECTRIC



SOFT TPE



ADJUSTABLE NOSEBRIDGE



Clear

Key Features

- ✓ Single wraparound polycarbonate lens
- ✓ Soft, ventilated TempFLEX™ temple tips provide a comfortable, secure fit
- ✓ Interchangeable and removable adjustable headband
- ✓ High performance TECTON 400 anti-scratch and anti-fog coatings applied to BOTH sides of the lens
- ✓ Soft, adjustable, anti-slip nose bridge
- ✓ Slim temples enhance comfort when worn with ear defenders
- ✓ Anti-reflective and high contrast lens options
- ✓ 100% metal-free
- ✓ Supplied with a microfibre storage pouch and premium cord



Blue

TECTON 400



Grey

TECTON 400



Amber

TECTON 400



LED

TECTON 400

 CAT II	Packaging	10 per box 200 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR	○	RLY00366	2C-1.2 e 1FTKN CE / Z87+ U6 IO 0	EN 166 FT CE / Z87+ AS/ NZS 1337.1 Lic, No. 1502	Impacts, UV radiations
GREY	●	RLY00367	5-2.5 e 1FTKN CE / Z87+ U6 L2.5 I 2	EN 166 FT CE / Z87+ AS/ NZS 1337.1 Lic, No. 1502	Impacts, UV and Solar radiations
BLUE	●	RLY00373	2-1.4 e 1FTKN CE / Z87+ U6	EN 166 FT CE / Z87+	Impacts, UV radiations
AMBER	●	RLY00368	2C-1.2 e 1FTKN CE / Z87+ U6	EN 166 FT CE / Z87+	Impacts, UV radiations
LED	●	RLY00372	2C-1.4 e 1FTKN CE / Z87+ U6 L1.3	EN 166 FT CE / Z87+ AS/ NZS 1337.1 Lic, No. 1502	Impacts, UV radiations

K&N Anti-scratch and
 Anti-fog / UV400 / Moisture repellent

• Product weight 23g | • Base lens 9mm



Riley®

Tilla™

Safety Glasses



LENGTH
ADJUSTABLE
TEMPLES



INCLINABLE
TEMPLES



OVER
SPECTACLES

Safety spectacle. Single piece polycarbonate frame suitable for over spectacle use.



Grey



Clear

Key Features

- ✓ Can be worn over spectacles
- ✓ Ventilated temple tips and side shields
- ✓ Brow protection
- ✓ Anti-fog lens
- ✓ UV400 protection
- ✓ Ideal to be worn by site visitors

CAT II CE	Packaging	10 per box / 200 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR ○	RLY00461	2C-1.2 e 1FT CE	e EN 166 FT CE / Z87+	Impacts, UV radiations	Anti-fog
GREY ●	RLY00462	5-2.5 e 1FT CE	e EN 166 FT CE / Z87+	Impacts, UV and Solar radiations	

• Product weight 42g | • Base lens 2mm



Riley® | Navigator™

Safety Glasses

Innovative yet iconic, the Riley Navigator offers a stylish solution to safety eyewear protection. The unique 4-base curve design combines classic contemporary sunglass styling with full safety features, allowing workers to experience retail design with protection in the workplace. Finished with a premium coating, the Navigator is available in a choice of 3 lens options including clear, polarised and grey.



Dark grey




Polarised
















Clear

Key Features

- ✓ Iconic fashion design
- ✓ 4 base curve lens
- ✓ Dual injected temples with soft tips
- ✓ Available in 3 lens options including clear grey and polarised
- ✓ Non-safety fashion option also available

	CAT II	Packaging	5 per box / 50 per case
		Certification	EN 166: 2001 PPE Regulation EU 2016/425

LENS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR	 RLY00526	2C -1.2 1F  / Z87+S U6 IO	 EN 166 F  / Z87+ AS/NZS 1337.1 LIC. SMK41066	Impacts, UV radiations	 K&N Anti-scratch and Anti-fog
DARK GREY	 RLY00522	5-3.1 1F KN  / Z87+S U6 L3 I3	 EN 166 F  / Z87+ AS/NZS 1337.1 LIC. SMK41066	Impacts, UV radiations	
POLARISED	 RLY00520	5-3.1 1F  / Z87+S U6 L3 I3	 EN 166 F  / Z87+ AS/NZS 1337.1 LIC. SMK41066	Impacts, UV radiations	Riley Anti-scratch and Anti-fog

• Product weight 30g | • Base lens 4mm



Riley Ready Reader™



Safety Glasses

The Riley READY READER is a sports style adjustable safety eyewear frame which is available in a range of magnifications from 1.00 to 3.00. The safety eyewear offers a comfortable and secure fit while providing additional magnification.



Key Features

- ✓ Available in 1.00, 1.50, 2.00, 2.50 and 3.00 magnifications
- ✓ Adjustable ratchet to position reading lens
- ✓ Dual-injected temples for enhanced comfort
- ✓ Anti-fog lens
- ✓ Soft flexible nose bridge for comfortable and secure fit
- ✓ Sports wraparound styling to compliment plano safety eyewear

CAT II CE	Packaging	10 per box / 200 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
RILEY READY READER 1.00	RLY00510	2C-1.2 e 1FT CE +1.0	e EN 166 FT CE	Impacts, UV radiations	Riley Anti-scratch and Anti-fog
RILEY READY READER 1.50	RLY00511	5-2.5 e 1FT CE +1.5	e EN 166 FT CE	Impacts, UV radiations	
RILEY READY READER 2.00	RLY00512	2-1.4 e 1FT CE +2.0	e EN 166 FT CE	Impacts, UV radiations	
RILEY READY READER 2.50	RLY00513	5-1.7 e 1FT CE +2.5	e EN 166 FT CE	Impacts, UV radiations	
RILEY READY READER 3.00	RLY00514	5-1.7 e 1FT CE +3.0	e EN 166 FT CE	Impacts, UV radiations	

• Product weight 44g



SPORT COLLECTION

Four sports-style prescription safety designs are available with four lens options: Polycarbonate, CR39, Trivex and Toughened Glass.

Riley® | Rokka™

Prescription strength up to a combined maximum of +/- 6.00



Frame Colour: Red ● Green ● Blue ● Grey ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
59mm	17mm	135mm	<input type="radio"/> Clear

Riley® | Script™

Prescription strength up to a combined maximum of +/- 4.00



Frame Colour: Black ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
60mm	18mm	130mm	<input type="radio"/> Clear

Riley® | Sivion™

Prescription strength up to a combined maximum of +/- 4.00



Frame Colour: Carbon Fibre Effect Black ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
61mm	18mm	125mm	<input type="radio"/> Clear

Riley® | Swarm™

Prescription strength up to a combined maximum of +/- 6.00



Frame Colour: Black/Blue ● Black/Grey ●

Black/Grey model features inner frame seal.

Lens (eyesize)	Nosepiece	Temple Length	Lens
59mm	18mm	130mm	<input type="radio"/> Clear



Classic COLLECTION

The Classic Collection includes eight contemporary prescription safety styles with slim, metal frames and transparent side shields.

Riley® | R101

Prescription strength up to a combined maximum of +/- 8.00



Frame Colour: Shiny Anthracite ●
Tip Colour: Solid Black ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
57mm	16mm	145mm	○ Clear

Riley® | R102

Prescription strength up to a combined maximum of +/- 8.00



Frame Colour: Brushed Anthracite ●
Tip Colour: Solid Black ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
54mm	15mm	140mm	○ Clear

Riley® | R103

Prescription strength up to a combined maximum of +/- 8.00



Frame Colour: Shiny Dark Navy ●
Tip Colour: Solid Black ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
55mm	16mm	140mm	○ Clear

Riley® | R104

Prescription strength up to a combined maximum of +/- 12.00



Frame Colour: Shiny Teal ● Gunmetal ● Navy (small) ●
Tip Colour: Solid Black ●

	Lens (eyesize)	Nosepiece	Temple Length	Lens
Small	48mm	15mm	130mm	○ Clear
Standard	52mm	15mm	130mm	○ Clear



Riley® | R105

Prescription strength up to a combined maximum of +/- 8.00



Frame Colour: Brushed Brown ●

Tip Colour: Solid Black ●

	Lens (eyesize)	Nosepiece	Temple Length	Lens
Standard	55mm	17mm	140mm	<input type="radio"/> Clear
Large	60mm	18mm	145mm	<input type="radio"/> Clear

Riley® | R106

Prescription strength up to a combined maximum of +/- 8.00



Frame Colour: Matt Black ●

Tip Colour: Solid Black ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
55mm	17mm	145mm	<input type="radio"/> Clear

Riley® | R107

Prescription strength up to a combined maximum of +/- 8.00



Frame Colour: Gunmetal ●

Tip Colour: Solid Black ●

Lens (eyesize)	Nosepiece	Temple Length	Lens
56mm	16mm	145mm	<input type="radio"/> Clear



Riley ^e | Reon™


Face Shield

The REON™ face shield provides the wear with 'Class A' face protection. High energy impact resistance to EN 166 'Level A' means that the face shield protects against impacts up to 190m/s (425mph). REON is also top of the class for comfort and vision clarity.



Key Features

- ✓ Fitted chinguard
- ✓ Ratchet adjustment allows for a secure fit to a range of head sizes
- ✓ Padded headband for a comfortable fit
- ✓ Adjustable pivot for wearer to select optimum angle
- ✓ Frameless design for a wide, unrestricted field of vision
- ✓ Simple, cost-efficient visor replacement system

 CAT II	Packaging	20 Polycarbonate per Carton 10 Acetate per Carton 20 Browguard per Carton
	Certification	EN 166: 2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
POLYCARBONATE AND CHINGUARD	RLY00351	2C-1.2 ^e 1AT 9 IO 0 CE	N/A	Impacts, UV radiations, molten metal splash	Anti-fog
BROWGUARD	RLY00352	N/A	^e EN166 AS/NZS 1337.1T CE		
ACETATE VISOR	RLY00355	2C-1.2 ^e 1FT3 CE +2.0	N/A	Impacts, UV radiations	

• Product weight 360g



Face Shield

The newly launched Alpha face shield protects against airborne droplets and offers direct splash protection. The Alpha face shield guarantees comfort thanks to the foam headband and adjustable elastic band. Lightweight and comfortable for all-day wear.



Key Features

- ✓ PETA visor
- ✓ Wide field of vision
- ✓ Can be worn over glasses
- ✓ Ultra-lightweight design
- ✓ Easy to fit by means of an elastic headband
- ✓ Metal-free
- ✓ Foam sweatband

	Item number	ASE9901
	Packaging	120 per case
	Certification	Tech spec based on EN 166: 2001 PPE Regulation EU 2016/425

ENHA 4430 Tomcat

Face Shield & Browguard

Face shield with headband. Headband (4430KD) can be adjusted with ratchet system.

Key Features

- ✓ Sweatband (443030), browguard (443010) + visor from Polycarbonate (4431), 1 mm or 1,5 mm, clear
- ✓ Suitable for: Laboratory, industry, electricians
- ✓ Material: PP: Browguard, polycarbonate: visor
- ✓ Colours: whale blue (browguard), visor clear (standard)
- ✓ High temp resistance: +55°C
- ✓ Low temp resistance: -5°C
- ✓ Weight: 320g / 350g

	Item number	4430PC000000
	Packaging	20 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425





Riley®

Velia™



OVER SPECTACLES

Safety Goggles

The Riley VELIA™ safety goggle has a clear polycarbonate wraparound lens offering 180° field of vision for complete eye protection.

Key Features

- ✓ Wide field of vision
- ✓ Padded face seal
- ✓ Strong polycarbonate lens
- ✓ Soft PVC/TPE frame and comfort fit nose bridge
- ✓ Quick release head strap
- ✓ Pivotal strap



Optional Easy-Change Buckle Headstrap



CAT II CE	Packaging	100 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR ○	RLY00221	2C-1.2 e 1BT9KN IMO 0 CE	e EN166 349BT AS/NZS CE	Impacts, Liquid Droplets, Molten Metal, Airborne Particles	TECTON 400™ K&N Anti-scratch and Anti-fog
BUCKLE HEADSTRAP	RLYA0004	N/A	N/A	N/A	N/A

• Product weight 84g

ENHA | EN 61 Goggle

Unvented Goggle

Key Features

- ✓ Unvented goggle, with adjustable headband
- ✓ Medium energy impact and liquid splash protection
- ✓ Material: Lens - polycarbonate, Frame - polycarbonate, woven elastic band
- ✓ Colours: Lens: clear, body; green, headband: black
- ✓ Storage: +5°C -40°C. Humidity less than 90%
- ✓ Anti-fog lens

CAT II CE	Packaging	200 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425



VERSIONS:	REF:	LENS MARKING:	FRAME MARKING:	PROTECTION AGAINST:	COATINGS:
CLEAR - AF ○	ENE00061RZ	2C - 1.2 ENHA 1BT CE / 787+	ENHA EN166 3 BT CE	Impacts, UV radiations	Anti-scratch and Anti-fog
CLEAR ○	ENE00061KZ	2C - 1.2 ENHA 1BT CE / 787+	ENHA EN166 3 BT CE	Impacts, UV radiations	Anti-scratch



Benchmark | BM30


Low Profile Panoramic Goggle

Key Features

- ✓ Excellent protection against liquid splash and dust
- ✓ Wide panoramic vision
- ✓ Universal nose bridge
- ✓ Soft TPR frame
- ✓ Indirect ventilation
- ✓ Adjustable head strap
- ✓ Optical class 1 Polycarbonate lens
- ✓ Can be worn over spectacles



Clear

CAT II 	Packaging	10 per box 100 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR ○	BME00030RZ	2C-1.2 ▶ 1BT 9 	▶ EN 166 3 4 9 BT  BM30	Impacts, UV radiations, dust, chemicals and molten metals	Anti-scratch and Anti-fog

Benchmark | BM08


Wraparound Style Safety Spectacle

Key Features

- ✓ Close-fitting protection
- ✓ Universal nose bridge
- ✓ Slim temple tips
- ✓ Anti-scratch lens
- ✓ Lightweight: only 28g
- ✓ Optical class 1 Polycarbonate lens



Clear

CAT II 	Packaging	200 per case
	Certification	EN 166: 2001 PPE Regulation EU 2016/425

VERSIONS	REF.	LENS MARKING	FRAME MARKING	PROTECTION AGAINST	COATINGS
CLEAR ○	RLYB0021	2C-1.2 ▶ 1FT 	▶ EN 166 FT  BM08	Impacts, UV radiations	Anti-scratch

• Product weight 23g | • Base lens 10mm



VIVID⁺ LENS CLEANING SOLUTIONS

Multi-Function Lens Cleaning Stations

Key Features

- ✓ Refillable dispenser unit
- ✓ Wall mountable
- ✓ Comprises 2 x 250 dry tissues and 500ml cleaning

Order Code: RLYA0018



Key Features

- ✓ Convenient 'ready-to-hang' solution
- ✓ Comprises 250 dry tissues and 500ml cleaning fluid
- ✓ Dispenser degradable after use

Order Code: RLYA0008



Dry Cleaning Tissue

Key Features

- ✓ Refillable dispenser unit
- ✓ Wall mountable
- ✓ Comprises 2 x 250 dry tissues and 500ml cleaning

Order Code: RLYA0007



Moist Lens Cleaning Tissues

Key Features

- ✓ Convenient 'ready-to-hang' dispensing solution
- ✓ Comprises 200 pre-moistened individually wrapped lens wipes
- ✓ Anti-static, anti-reflective and anti-bacterial
- ✓ Alcohol and silicone-free
- ✓ Dispenser degradable after use

Order Code: RLYA0005



Lens Cleaning Spray

Key Features

- ✓ 500ml spray bottle
- ✓ Alcohol and silicone-free
- ✓ Anti-reflective
- ✓ Anti-static

Order Code: RLYA0006



Anti-Fog Wipes

Key Features

- ✓ Safely and quickly cleans and prevents fogging on eyewear lenses, electronic surfaces and other optics
- ✓ Convenient 'ready-to-hang' dispensing solution
- ✓ Comprises 200 pre-moistened individually wrapped lens wipes
- ✓ Anti-static, anti-reflective and anti-bacterial
- ✓ Alcohol and silicone-free
- ✓ Dispenser degradable after use

Order Code: RLYA0027





Contact our UK team

Need more information?

Discuss your requirements with one of our
healthcare PPE specialists.

Our UK team is here to help and advise you.

www.globusgroup.com

T: +44 (0)161 877 4747

E: healthcare@globusgroup.com



Hand Protection



Respiratory Protection



Head Protection



Hearing Protection



Eye Protection



Protective Clothing